






CATALOGUE

Ultra-Clean, Reliable Sealing for Analytical Laboratories



NOKE
LAB

Your Trusted Manufacturer for High-Quality Chromatography Consumables

MICRO SIZE, PURE KEEPS

Vials and Inserts

Closures and Septa,

Filters and SPE Columns

www.nokesci.com



MICRO SIZE, PURE KEEPS!

Your Trusted Manufacturer for
High-Quality Chromatography Consumables

Vials and Inserts, Closures and Septa, Filters and SPE Columns




www.nokesci.com

NOKE LAB - Micro Size, Pure Keeps

Tiny Parts. Mighty Manufacturing. Pure Chromatography Beyond.

At **NOKE LAB**, we specialize in engineering precision chromatography consumables that deliver uncompromised analytical performance.



From **GC/HPLC vials, septa, and closures** to **syringe filters**, our tiny yet vital components are built with uncompromising manufacturing excellence to ensure complete contamination control—even in your most demanding applications.

Manufacturing Excellence. Trusted Quality. Proven Performance.

Manufacturing Excellence

10,000 sqm state-of-the-art facility featuring:
5,500 sqm precision assembly workshop
5,500 sqm injection molding center
Class 100,000 cleanroom environments

Fully automated production lines ensuring:

Batch-to-batch consistency
Human contamination elimination
24/7 manufacturing capability

The NOKE LAB Difference

- Ultra-Clean Production - Particle-free consumables for sensitive MS detection
- Robotic Precision - Dimensional accuracy for modern autosamplers
- Rapid Customization - 24-hour prototype turnaround
- Global Logistics - Reliable just-in-time delivery
- Innovation Partnerships - Joint R&D with top technical universities

Applications

Regulated pharmaceutical analysis

Trace environmental testing

Food safety monitoring

Advanced research methodologies

Where micro-scale precision meets macro-level reliability

Global Presence & Market Recognition

NOKE LAB serves over 50 countries across Europe, North America, South Asia, the Middle East, and beyond. Trusted by laboratories worldwide, NOKE has become a recognized name in chromatography consumables, known for quality, precision, and reliability.

We are committed to enabling scientific breakthroughs and supporting global diagnostics through Continuous innovation | Strict quality control | Reliable service and support
NOKE LAB is proud to play a part in discoveries that shape a healthier, more advanced world.



Bonded Closures and Certified Vial Kits

- UltraBond Closures ND9 /1
- UltraBond Closures ND24 /1
- LC/MS and GC/MS Certified Vial Kits /2

Screw Neck ND8

- PP Screw Closures ND8 with Septa /3
- Screw Neck Vials ND8 and Micro Inserts for ND8 Vials /3

Short Thread Screw Neck ND9

- PP Screw Closures ND9 /4-5
- UltraBond Closures ND9 /5
- Short Thread Screw Neck Vials ND9 and Micro Inserts for ND9 Vials /5-6
- Micro-Vials ND9 for Small Sample Volumes (High Recovery Vials) /6
- Pre-Sealed Vials with Closures ND9 /7
- LC/MS and GC/MS Certified Vial Kits /8

Screw Neck ND10

- PP Screw Closures ND10 with Septa /8
- Short Thread Screw Neck Vials ND10 /9

Crimp Neck ND11

- Sliver Crimp-top Aluminum Cap /9
- Crimp Neck Vials ND11 /9
- Micro-Vials ND11 for Small Sample Volumes (High Recovery Vials) /10

Snap Ring ND11

- PP Snap Closures ND11 with Septa /10
- Snap Ring Vials ND11 /11

Screw Neck ND13

- PP Screw Closures ND13 with Septa /11
- Screw Neck Vials ND13 /11

Screw Neck ND15

- PP Screw Closures ND15 with Septa /12
- Screw Neck Vials ND15 /12

Headspace Screw Neck ND18

- Magnetic Screw Closures ND18 with Septa /12
- Precision Thread Screw Neck Vials ND18 /13

Headspace Crimp Neck ND20

- Aluminum Crimp Closures and Magnetic Crimp Closures ND20 with Septa /13-14
- Crimp Neck Vials ND20 /14-15

Screw Neck ND24 (EPA)

- PP Screw Closures ND24 /16
- UltraBond Closures ND24 /17
- Screw Neck Vials ND24 for TOC / VOC Analysis and Sample Storage /17
- TOC Sets: Pre-Cleaned and Pre-Sealed Vials with Closures ND24 /18
- Per-Seals Kits ND24 /18

Diverse Vials / Tubes / Bottles and Caps

- COD Vials and Closures / Pre-Sealed Kits /19
- Shell Vials and Closures /19
- Sample Cell Vials with Closures /20
- Scintillation Vials with Closures /20
- NMR Tubes /20
- Blue Screw Cap Reagent Bottles (GL45 / GL80) /21
- HPLC Safety Caps & Waste Container Caps /21-22

Tools and Racks for Chromatography

- Racks for Vials /23
- Manual Crimpers and Decappers (11 mm / 20 mm) /23
- Electronic Crimpers and Decappers /23

Syringe Filters and Membranes

- Nylon (NY) /25
- PTFE (Hydrophilic/Hydrophobic) /26
- PVDF (Hydrophilic/Hydrophobic) /26
- Polyethersulfone (PES) /25
- Mixed Cellulose Esters (MCE) /25
- Sterilized Syringe Filters /28
- Syringless Filters Vial /29-30
- Membranes (25mm / 47 mm / 50 mm) /31-32

SPE Cartridges

- Polymer-Based Sorbents SPE /33-35
- Silica-Based Sorbents SPE /33-35
- Adsorption Sorbents SPE /33-35
- Mixed-Mode Sorbents SPE /33-35

General Use Plastic Consumables

- Liquid Handling and Storage / Molecular Biology /37-46
 - Cell Culture Series Products /47-50
 - Others
-

Chromatography Vial

As a trusted partner in the chromatography field, NOKE LAB has been designing and manufacturing vials, caps, and seals since the early days of automated analytical techniques.

Our deep expertise enables us to deliver optimized consumable combinations fully compatible with today's most widely used instruments.



- Pass the ISO9001: 2015 quality certification. Audited and issued by the most stringent SGS in the industry.
- Experienced technician with equipped R&D center to fulfill new research and product development.
- Be 100% compatible with original Agilent, Waters, Virian, Shimadzu, etc.
- All materials are produced in a class 100,000 cleanroom.
- OEM services are available.

Proper septa selection is critical to preventing contamination and avoiding damage to the needle and instrument.

1.5ml

1.8ml

2ml



1.5mL Vials
Usable Volume

What is glass 7.0/5.0/3.3 type?

Type	Material	Coefficient of Thermal Expansion	Features
7.0	Borosilicate Glass	$(6.2\sim 7.5) \times 10^{-6} K^{-1} (20\sim 300\text{ }^{\circ}C)$	Economy Type Competitive Price For Disposable Quick Test
5.0	Borosilicate Glass	$\leq 5 \times 10^{-6} K^{-1} (20\sim 300\text{ }^{\circ}C)$	higher quality Type For Test and Save Samples
3.3	Borosilicate Glass	$3.3 \times 10^{-6} K^{-1} (40\sim 500\text{ }^{\circ}C)$	Resistant to strong acids and alkalis, with no metal ion leaching.

Guide to Cap Selection

Choosing the right vial cap is key to protecting sample integrity, avoiding contamination, and matching chromatographic needs. Below is a concise guide by cap type and use case:

Screw Caps

Use For: LC/LC-MS, routine analysis, high-throughput auto testing (works with auto-samplers from Agilent, Shimadzu, Waters, etc.), frequent opening/closing.

Advantages: Tool free (twist by hand/gloves), good resealability, fits standard threads (8-425, 13-425).

Recommended Gaskets: PTFE/silicone (versatile for most solvents); pre-slit PTFE/silicone (eases auto-sampler needle use, reduces vacuum).



Crimp Caps

Use For: GC/GC-MS, high-security needs (food testing, forensics—prevents tampering), volatile organic compound (VOC) analysis.

Advantages: Excellent sealing (stays airtight during temperature changes), good for long storage/high-temperature tests.

Recommended Gaskets: PTFE/butyl rubber (chemical-resistant, withstands repeated punctures); PTFE/silicone /PTFE (no silicone shedding, ideal for trace analysis).



Snap Caps

Use For: Short-term storage (non-volatile samples), low-frequency use (e.g., small-lab manual sampling).

Advantages: Tool-free, quick to assemble/disassemble, cost-effective for basics.

Note: Not for high-precision tests or volatile sample storage (weaker sealing).



Quick Selection Tips

Volatile compounds/GC/GC-MS: Choose crimp caps + PTFE-lined gaskets (prevents sample loss).

High-throughput LC/LC-MS: Pick screw caps + pre-slit gaskets (boosts auto-sampler efficiency).

Trace analysis/food/forensics: Use crimp caps + anti-shedding gaskets (e.g., PTFE/silicone/PTFE).

Corrosive solvents: Match caps with pure PTFE gaskets (enhances chemical resistance).

SEPTA COMPATIBILITY

Reagent	PTFE	PTFE/ Silicone rubber	PTFE/Silicone rubber/ PTEE	PTEE/ Red rubber	Fluorine rubber	PTFE/Butyl rubber
Acetonitrile	●	●	●	●		●
Hexane, Heptane, Methane	●		●	●	●	
Methanol	●	●	●	●		●
Benzene	●		●		●	
Tetrahydrofuran	●		●			
Toluene	●		●			
Dimethylformamide DMF	●	●	●			●
Dimethyl Sulfone DMSO	●	●	●			●
Ether	●	●	●			
Chloromethane Chloride	●		●		●	
Ethanol	●	●	●	●	●	●
Acetic Acid	●	●	●			●
Acetone	●	●	●			
Phenol	●	●	●		●	●
cyclohexane	●		●	●	●	

VIALS COMPATIBILITY

Manufacturer	Autosampler	9mm Screw	13-425 Screw	11mm Snap&Crimp	HeadSpace	ND24 Screw
Agilent	1100/1200/1260Infinity/1290 Infinity	●		●		
	G1888A / HS7694 / HS7694		●		●	
	Agilent 7673A	●		●		
	Agilent 7683A	●		●		
	7695A / Tekmar SOLATEk72 /Archon Pluge + Trap / AQUATEk 70					●
	79855(A) / 5880 / 5890 /6850 (27 Pos. Tray) / 6890	●		●		
	CTC HTS + HTC PAL + CTC GC PAL	●		●		
	CTC Combi PAL				●	
	7693A		●	●		
	7697A				●	
	7650A	●		●		
Waters	ACQUITY UPLC Systeme	●	●			
	Wisp 48 position		●			
	Wisp 96 position / 717, 96 Position Carousel					
	717 plus	●	●	●		
	717, 48 Position Carousel		●			
	Alliance	●		●		
	Alliance GPC 2000		●		●	
	Alliance HT Syst. / Alliance 2790 / 2795 /Alliance 2690 / 2695	●		●		
	Acquity Sample Organizer / Acquity / CapLC /Waters / Nano Acquity	●				
	Acquity H-Class / Alliance HTS					
	Model 2767 / Model 2707 / Model 2777	●		●		
Shimadzu	AOC-5000	●		●	●	
	AOC-14/1400/AOC-17/AOC-20/20i/20s 150 Pos. Tray	●	●	●	●	
	AOC-20/20i/20s 96 Pos. Tray		●			
	LC-20A/SIL-6B/SIL-7A/SIL-8A/SIL-9A	●	●	●		
	SIL-2AS/SIL-6A	●	●	●		
	SIL-10A/SIL-10AF/SIL-10AP/ SIL-10Ai /SIL-10AxL	●		●		
	SIL-10HTA/SIL-10HTC 140 Pos. Tray /LC2010C+LC2010A 140 Pos. Tray	●		●		
	SIL-10ADvp	●	●	●		
	HTA 200 H/ HSS-2B				●	
	LC2010C + LC2010A 100 Pos. Tray		●			
	ASI-V					●
	SIL 30-ACMP	●		●		
	SIL-20AxR/ SIL-20ACxR (Prominence)	●	●	●		
	SIL-20AHT, SIL-20AHT	●	●	●		
PerkinElmer	Series 200, 85 vial tray			●		
	Series 200, 81/ 100 vial tray			●		
	Series 200, 205 vial tray/ ISS-225, 205 vial tray			●		
	AI-1/ AS-100/ AS-100B			●		
	AS2000/ AS2000B			●		
	AS-300/ AS8300/ Autosystem			●		
	TurboMatrix HS16/HS40/ HS40 xL/HS40Trap/HS110/HS110 Trap			●		
	ISS-100, 85 vial tray/ISS-200, 85 vial tray/ISS-225, 100 vial tray+80 vial tray			●		
	Clarus 400, 500, 600/590 GC, 690 GC			●		
	Flexar Fx-20 UHPLC; LC	●		●		

VIALS COMPATIBILITY

Manufacturer	Autosampler	9mm Screw	13-425 Screw	11mm Snap&Crimp	HeadSpace	ND24 Screw	
CTC (LEAP)	LC PAL (216 Pos.)	●		●	●		
	HTx PAL/Combi PAL						
	Combi PAL (98 Pos. Tray), GC PAL (98 Pos. Tray)	●		●			
	Combi PAL SPME Mode (98 Pos. Tray)	●		●			
CTC	PAL HPLC-Systems/PAL Combi-xt Liqid Mode/PAL HTC-xt	●		●	●		
	Combi-xt Headspace Option			●	●		
	GC-xt Headspace Option/A200 LC	●		●	●		
	Combi-xt SPME Options			●	●		
	A200S	●		●			
	HS 500				●		
	PAL RTC/ PAL RSI				●		
	PAL LSI	●	●	●			
	Thermo Fisher	UltiMate Analytical, conical, WPS-3000 SL			●	●	
		UltiMate™ Micro, conical, WPS-3000 SL				●	
UltiMate™ Semipreparative, WPS-3000 SL			●		●		
AS-AP* (120 Pos. 1.5 mL), (3 x plates)		●		●			
AS2000 30 vial tray/ HS250/ HS500/HS800/ HS2000					●		
AS2000 90 vial tray (Trace GC)		●		●			
AI3000 (II)/ AS3000 (II)/ AS3500(Trace GC + Focus GC)		●		●	●		
A200LC/ Accela™ OpenAutosampler		●		●			
A200S		●		●			
AS800, 42 vial tray		●		●			
AS800, 60 vial tray		●		●			
Dionex™ AS-AP		●		●			
Dionex™ UltiMate WPS-3000		●	●	●			
Dionex™ AS 40			●				
Trace™ 1300 Series/ Trace™ 1310 Series /AI/AS 1310Series /TriPlus™ 100 LS		●	●	●			
UltiMate™ WPS-3000RS		●		●	●		
Dionex™ AS-AP		●		●			
Dionex™ UltiMate WPS-3000		●	●	●			
Dionex™ AS 40			●				
VWR		Hitachi Chromaster	●	●	●		
	HPLC-System Primaide	●	●	●			
	L2200/L2200-U/L7200/L7250/655-A40/L-9100	●					
	L2200 (LaChrom Elite) (128 Pos. Tray)		●				
	5210 (Chromaster)	●	●	●			
Gilson	201/202/221/222/Aspec	●	●				
	231/401/232/402/Aspec xli/ Aspec xL4	●					
	Nano Injektor	●					
	235/235P/SP 235/SP 235P	●					
Beckman	501/502/502e/507/ 507e	●		●			
	508 (System Gold)	●			●		
	Marathon/Promis	●		●			
	Triathlon, Standard Tray	●		●	●		
	Triathlon, LSV Tray		●				
	Triathlon, Super-LSV Tray				●		

Autosampler compatibility and cross-references

Series	Noke	Agilent	MN	Waters	Thermo	Shimadzu	PE
ND8 Screw Vials	NKSV-2-8	5183-4428					
	NKSV-2-8-A	5183-4429					
	NKSV-2-8-W				60180-508	038-00165-06	
	NKSV-2-8-A-W				60180-560		
ND8 Caps & Septum	NKSSC-8-B	5183-4442	70245		60180-665		
ND8 Inserts	NKI-250	5183-2090	702005				
	NKCI-150	5183-2089	702968.1				
	NKIF-150	5183-2088	70284				
ND9 Screw Vials	NKSV-2-9	5182-0714	702282	186000273	60180-723	980-08337	N9306201
	NKSV-2-9-W	5182-0715	702283	186000273	60180-509	980-08337	N9306201
		5183-2068					
	NKSV-2-9-A-W	5182-0716	702284	186000848	60180-561		N9306220
	NKSV-2-9-A-U	5188-6535	702293		60180-724		
		5188-6536		186002626			
	NKPSV-03-9-C	5190-2242	702135				
		5190-2243					
	NKSV-03-9	5188-6591	702007				
NKSV-03-9-A	5188-6592	702008					
NKSV-13-9	5183-2030	702006		186002802			
NKSV-13-9-A	5183-2073						
ND9 Screw Vials	NKSSC-9-A	5182-0717			60180-728		
	NKSSC-9-B	5182-0720	702287.1		60180-516		N9306202
	NKSSC-9-F-P	5183-2076			60180-672		
	NKSSC-9-B2	5185-5823		186000274			
	NKSSC-9-C	5182-0723			60180-670		
	NKSSC-9-B2-P	5185-5824		186000305	60180-674		
ND9 Inserts	NKI-300	5181-3377	702825				
	NKCI-250-I	5183-2085	702716				
	NKIF-250-I	5181-1270	702818	WAT094170			
ND11 Crimp Vials	NKCV-2-11-W-U	5182-0543	702885	WAT094222	60180-502		N9301385
	NKCV-2-11-A-W-U	5181-3376	702892	WAT094223			N9302680
	NKCV-2-11-U	5181-3375	70201HP				
ND11 Crimp caps & Septum	NKSSCC-11-B	5182-0552	70288	PSL404219	60180-706		N9306228
	NKSSCC-M-11-R	5181-1210			60180-705		N9306015
ND11 Snap Vials	NKNV-2-11-W	5182-0546	702713	WAT094219	60180-710		N9303418
	NKNV-2-11-A-W	5182-0545	702712	WAT094220			N9306208
ND11 Snap caps & Septum	NKSSNC-11-B	5182-0541	702710	186002649	60180-713		N9303419
	NKSSNC-11-A	5182-0564			60180-676		
1ML Shell Vial	NKSV-1			WAT025054C			
ND13 Screw Vials	NKSV-4-13	5183-4448	702962	186000840		220-90393-00	N9306247
	NKSV-4-13-A	5183-4450	702973				
ND13 Caps & Septum	NKSSC-13-B	5183-4464		186000841		220-91521-14	
	NKSC-13	5183-4461		WAT072711		220-90394-00	
ND18 Headspace Vials	NKSV-10R	5188-5392	702866			220-94796-07	N6356479
	NKSV-10R-A	5188-6538					
	NKSV-20R	5188-2753	702826			220-94796-01	N9306075
	NKSV-20R-A	5188-6537	702826.2				
ND18 Caps & Septum	NKSSC-18M-F	5188-2759	702055				
	NKSSC-18M-T		702827				
ND20 Headspace Vials	NKCV-10F	5182-0838	702918				
	NKCV-20F	5182-0837	702261				
	NKCV-20R	5183-4474	702263		60180-506		N9306079
	NKCV-20R-A	5190-2239	702263.2				
ND20 Crimp caps & Septum	NKCC-20	9301-0721			60180-512		
	NKCC-20S	9301-0718					
	NKSS-GB-20	5183-4476					
	NKSSCC-20-G	5183-4477			60180-511		
	NKSSCC-20S-G	5183-4478	702835				
	NKSSCC-20M-G		702837		60180-520		
ND24 Screw Vials	NKSV-40-E		702023				
	NKSV-40-E-A		702024				
ND24 Caps & Septum	NKSSC-22-G	5183-4308	702058			038-00165-40	N6352033

Bonded Closures

Bonded cap and septa are produced without adhesives or chemicals, using advanced bonding machines to effectively eliminate contamination risks.

Molecular-Level Bonding · Zero Leakage Risk

Adopting high-temperature and high-pressure hot-press bonding technology, it achieves seamless integration of the PTFE surface layer, silicone functional layer, and PP support layer, completely eliminating the problems of "delamination and peeling" that plague conventional gaskets.



Exceptional Chemical Inertness · Protect Sample Purity

Its surface uses high-purity PTFE ($\geq 99.95\%$ purity) that resists strong acids, alkalis and organic solvents (methanol, acetonitrile, n-hexane, concentrated HCl, etc.), with no extractables to taint samples. Paired with the elastic seal of silicone/butyl rubber layer, it prevents volatilization and external impurities, ensuring no interfering peaks in LC/MS trace analysis.

UltraBond Closures ND9



NKSSC-9-A2



NKSSC-9-A2-P



NKSSC-9-B2



NKSSC-9-B2-P

Description	Septum	Slit Type	Cat. No.	Pack
Bonded, 9mm Blue Short Screw Polypropylene Cap with Septa, 5.5mm Centre Hole	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKSSC-9-A2	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKSSC-9-A2-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKSSC-9-B2	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKSSC-9-B2-P	100pcs/bag

ND24 Bonded Caps with Septas



NKSSC-22B-H



NKSSC-22B-S



NKSSBC-22B-H



NKSSBC-22B-S



NKSSBC-22B-S-15

Description	Septum	Cat. No.	Pack
Bonded, 24mm White Screw Polypropylene Cap with Septa, 15mm Centre Hole	Nature PTFE/ Nature Silicone, 22*3.0mm	NKSSC-22B-H	100pcs/bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22B-S	100pcs/bag
Bonded, 24mm White Screw Polypropylene Cap with Septa, 16mm Centre Hole	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSBC-22B-H	100pcs/bag
	Nature PTFE/ white Silicone, 22*3.0mm	NKSSBC-22B-S	100pcs/bag
	Nature PTFE/ white Silicone, 22*1.5mm	NKSSBC-22B-S-15	100pcs/bag

Certified Vial Kits

LC/MS and GC/MS Certified Vial Kits

Certification Support: Every Piece Meets Professional Standards

Full batches have passed USP Class VI Biocompatibility Certification, ISO 9001 Quality Management System Certification and RoHS Environmental Certification. Complete batch test reports (including extractable testing and particle content analysis) are provided.

Among them, USP Class VI Certification ensures no harmful substance migration when in contact with food and pharmaceuticals, fully meeting the sample processing and storage requirements under GLP standards.



NKCCSV-2-9-A



NKCCSV-2-9-B



NKCASV-2-9-A



NKCASV-2-9-B



NKCCSV-2-9-A-P



NKCCSV-2-9-B-P



NKCASV-2-9-A-P



NKCASV-2-9-B-P

Description	Septum	Slit Type	Cat. No.	Pack
Certified Kits, Pre-sealed, 2mL Amber Screw Vial w/Marking Spot+Blue Cap	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKCASV-2-9-A	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKCASV-2-9-A-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKCASV-2-9-B	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKCASV-2-9-B-P	100pcs/bag
Certified Kits, Pre-sealed, 2mL Clear Screw Vial w/Marking Spot+Blue Cap	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKCCSV-2-9-A	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKCCSV-2-9-A-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKCCSV-2-9-B	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKCCSV-2-9-B-P	100pcs/bag

8mm Black Short Screw Polypropylene Cap

Compatible with autosamplers of major brands such as Shimadzu, Agilent, Waters, and Varian. Widely applied in fields including pharmaceutical R&D and quality control, environmental pollutant detection, food additive analysis, and chemical raw material purity testing. It meets the requirements for sample storage, pretreatment, and injection in chromatographic analyses such as HPLC and GC.



Description	Septum	Slit Type	Cat. No.	Pack
8mm Black Screw Polypropylene Cap, 5.5mm Centre Hole	White PTFE/Red Silicone	non-slit	NKSSC-8-A	100pcs/bag
	White PTFE/Red Silicone	"-" Pre-slit	NKSSC-8-A-P	100pcs/bag
	Red PTFE/White Silicone	non-slit	NKSSC-8-B	100pcs/bag
	Red PTFE/White Silicone	"-" Pre-slit	NKSSC-8-B-P	100pcs/bag
	Blue PTFE/White Silicone	non-slit	NKSSC-8-F	100pcs/bag
	Blue PTFE/White Silicone	"-" Pre-slit	NKSSC-8-F-P	100pcs/bag

Screw Neck Vials ND8 and Micro Inserts for ND8 Vials



Description	Material	Writing area	Cat. No.	Pack
2mL Clear Screw Neck Vial, 8-425 Thread, 11.6*32mm	Boro 7.0	NO	NKSU-2-8	100pcs/box
	Boro 7.0	YES	NKSV-2-8-W	100pcs/box
	Boro 5.0	NO	NKSU-2-8-U	100pcs/box
	Boro 5.0	YES	NKSV-2-8-W-U	100pcs/box
2mL Amber Screw Neck Vial, 8-425 Thread, 11.6*32mm	Boro 7.0	NO	NKSV-2-8-A	100pcs/box
	Boro 7.0	YES	NKSV-2-8-A-W	100pcs/box
	Boro 5.0	NO	NKSV-2-8-A-U	100pcs/box
	Boro 5.0	YES	NKSV-2-8-A-W-U	100pcs/box
Glass Micro-Insert for ND8 Vials suits for 8-425 Vials	150uL Mandrel Interior & Polymer Feet, 29*5mm		NKIF-150	100pcs/bag
	150uL Conical Micro-Insert, 31*5mm		NKCI-150	100pcs/bag
	250uL Micro-Insert, 31*5mm, Flat Bottom		NKI-250	100pcs/bag

PP Screw Closures ND9

- The ND9 vials can be used on all common autosamplers that perfectly compatible with Agilent, Waters, Thermo Scientific, etc.
- Wide opening improves sample accessibility and reduces the risk of autosampler needle damage.
- Optional white writing patch and fill marks for easy labeling.
- Made from borosilicate glass to meet different customer requirement.
- Pre slit shape are -slit, Y-shape, cross-shape available.



NKSSC-9-A



NKSSC-9-A-CP



NKSSC-9-A-P



NKSSC-9-A-YP



NKSSC-9-B



NKSSC-9-B-CP



NKSSC-9-B-P



NKSSC-9-B-YP



NKSSC-9-F



NKSSC-9-F-CP



NKSSC-9-F-P



NKSSC-9-F-YP

Description	Septum	Slit Type	Cat. No.	Pack
9mm Blue Short Screw Polypropylene Cap with Septa, 5.5mm Centre Hole	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKSSC-9-A	100pcs/bag
		Cross Pre-slit	NKSSC-9-A-CP	100pcs/bag
		"-" Pre-slit	NKSSC-9-A-P	100pcs/bag
		"Y" Pre-slit	NKSSC-9-A-YP	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKSSC-9-B	100pcs/bag
		Cross Pre-slit	NKSSC-9-B-CP	100pcs/bag
		"-" Pre-slit	NKSSC-9-B-P	100pcs/bag
		"Y" Pre-slit	NKSSC-9-B-YP	100pcs/bag
	Blue PTFE/White Silicone, 9*1.0mm	non-slit	NKSSC-9-F	100pcs/bag
		Cross Pre-slit	NKSSC-9-F-CP	100pcs/bag
		"-" Pre-slit	NKSSC-9-F-P	100pcs/bag
		"Y" Pre-slit	NKSSC-9-F-YP	100pcs/bag

PP Magnetic Screw Closures ND9/ UltraBond Closures ND9



NKSSC-9M-A-1



NKSSC-9M-B-1



NKSSC-9-A2



NKSSC-9-A2-P



NKSSC-9-B2



NKSSC-9-B2-P

PP Screw Closures ND9/ UltraBond Closures ND9

Description	Septum	Slit Type	Cat. No.	Pack
ND9 9mm golden Magnetic Open top Screw Blue cap with Septa, Central hole 6.0mm	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKSSC-9M-A-1	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKSSC-9M-B-1	100pcs/bag
Bonded, 9mm Blue Short Screw Polypropylene Cap with Septa, 5.5mm Centre Hole	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKSSC-9-A2	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKSSC-9-A2-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKSSC-9-B2	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKSSC-9-B2-P	100pcs/bag

Short Thread Screw Neck Vials ND9 and Micro Inserts for ND9 Vials

Compatible with all common autosamplers. They are perfectly match used on Agilent, HTA, Shimadzu, Thermo, Waters® etc instruments.

The mainstream neck size for 2mL vials, paired with 9mm bonded septa (e.g., 9-425). Features a threaded seal design, enabling fast and precise puncture by autosamplers.



NKSV-2-9
NKSV-2-9-U
NKSV-2-9-U1

NKSV-2-9-W
NKSV-2-9-A-W-U
NKSV-2-9-W-U1

NKSV-2-9-A
NKSV-2-9-A-U
NKSV-2-9-A-U1

NKSV-2-9-A-W
NKSV-2-9-A-W-U
NKSV-2-9-A-W-U1



Description	Material	Writing area	Cat. No.	Pack
2mL Clear Short Thread Vial, 9mm Thread, 11.6*32mm	Boro 7.0	NO	NKSV-2-9	100pcs/box
	Boro 7.0	YES	NKSV-2-9-W	100pcs/box
	Boro 5.0	NO	NKSV-2-9-U	100pcs/box
	Boro 5.0	YES	NKSV-2-9-W-U	100pcs/box
2mL Amber Short Thread Vial, 9mm Thread, 11.6*32mm	Boro 7.0	NO	NKSV-2-9-A	100pcs/box
	Boro 7.0	YES	NKSV-2-9-A-W	100pcs/box
	Boro 5.0	NO	NKSV-2-9-A-U	100pcs/box
	Boro 5.0	YES	NKSV-2-9-A-W-U	100pcs/box
2mL Silanized Clear Short Thread Vial , 9mm Thread, 11.6*32mm	Boro 5.0	NO	NKSV-2-9-U1	100pcs/box
	Boro 5.0	YES	NKSV-2-9-W-U1	100pcs/box
2mL Silanized Amber Short Thread Vial , 9mm Thread, 11.6*32mm	Boro 5.0	NO	NKSV-2-9-A-U1	100pcs/box
	Boro 5.0	YES	NKSV-2-9-A-W-U1	100pcs/box

Short Thread Screw Neck Vials ND9 and Micro Inserts for ND9 Vials



Description	Material	Writing area	Cat. No.	Pack
Clear PP Short Screw Micro-Vial, 9mm Thread, 11.6*32mm	Volume : 2 ml		NKPSV-2-9	100pcs/box
	Volume : 0.3 ml		NKPSV-03-9-C	100pcs/box
Amber PP Short Screw Micro-Vial, 9mm Thread, 11.6*32mm	Volume : 2 ml		NKPSV-2-9-A	100pcs/box
	Volume : 0.3 ml		NKPSV-03-9-A	100pcs/box
Glass Micro-Insert for 9-425 vials	Shape: 250ul Conical Bottom, 31*6mm (box)		NKCI-250-I	100pcs/box
	Shape: 250ul Conical Bottom, 31*6mm (bag)		NKCI-250-I	100pcs/bag
	Shape: 300ul Flat Bottom, 31*6mm (box)		NKI-300	100pcs/box
	Shape: 300ul Flat Bottom, 31*6mm(bag)		NKI-300	100pcs/bag
	Shape: 250uL Mandrel Interior & Polymer Feet, 29*5.7mm (box)		NKIF-250-I	100pcs/box
	Shape: 250uL Mandrel Interior & Polymer Feet, 29*5.7mm(bag)		NKIF-250-I	100pcs/bag

Micro-Vials ND9 for Small Sample Volumes (High Recovery Vials)

Noke high recovery vials are the vial of choice when working with limited sample volumes but not wanting to manually combine a standard 2 mL vial with a separate vial insert. They ensure you extract the maximum possible amount of your precious sample from the autosampler needle.

Description		Volume	Cat. No.	Pack
ND9 Clear Vial High recovery		0.3ml Integrated w/Micro-Insert Base Bonded	NKSV-03-9	100pcs/box
ND9 Amber Vial High recovery			NKSV-03-9-A	100pcs/box
ND9 Clear Vial High recovery		1.1ml vial, 25µL Reservoir	NKSV-11-9	100pcs/box
ND9 Amber Vial High recovery			NKSV-11-9-A	100pcs/box
ND9 Clear Vial High recovery		1.3ml, 30µL Reservoir	NKSV-13-9	100pcs/box
ND9 Amber Vial High recovery			NKSV-13-9-A	100pcs/box
ND9 Clear Vial High recovery		1.5ml champagne Shape base	NKSV-15-9	100pcs/box
ND9 Amber Vial High recovery			NKSV-15-9-A	100pcs/box

Pre-Sealed Vials with Closures ND9

- Transparent short-thread glass vial: 2ml capacity, 9-425 specification, made of high-transparency material for clear observation of sample status, with printed markings for easy traceability;
- Blue open-top cap + composite septum: PTFE/silicone dual-material composite (9*1.0mm, no pre-slit), featuring secure sealing and contamination resistance;
- Non-screw cap design: Compatible with mainstream chromatographic instruments, enabling convenient operation and saving experimental time.



Description	Vials	Septum	Slit Type	Cat. No.	Pack
ND9 Sampler Vials Kit (unassembled) Clear vial +blue cap	2ml Clear vials w/writing spot	White PTFE/Red Silicone	non-slit	NKCSVK-2-9-A	100pcs/box
			"-" Pre-slit	NKCSVK-2-9-A-P	100pcs/box
		Red PTFE/White Silicone	non-slit	NKCSVK-2-9-B	100pcs/box
			"-" Pre-slit	NKCSVK-2-9-B-P	100pcs/box
		Blue PTFE/White Silicone	non-slit	NKCSVK-2-9-F	100pcs/box
			"-" Pre-slit	NKCSVK-2-9-F-P	100pcs/box
ND9 Sampler Vials Kit (unassembled) Amber vial +blue cap	2ml Amber vials w/writing spot	White PTFE/Red Silicone	non-slit	NKASVK-2-9-A	100pcs/box
			"-" Pre-slit	NKASVK-2-9-A-P	100pcs/box
		Red PTFE/White Silicone	non-slit	NKASVK-2-9-B	100pcs/box
			"-" Pre-slit	NKASVK-2-9-B-P	100pcs/box
		Blue PTFE/White Silicone	non-slit	NKASVK-2-9-F	100pcs/box
			"-" Pre-slit	NKASVK-2-9-F-P	100pcs/box

Description	Vials	Septum	Slit Type	Cat. No.	Pack
ND9 Sampler Vials Kit (assembled) Clear vial +blue cap	2ml Clear vials w/writing spot	White PTFE/Red Silicone	non-slit	NKCSV-2-9-A	100pcs/box
			"-" Pre-slit	NKCSV-2-9-A-P	100pcs/box
		Red PTFE/White Silicone	non-slit	NKCSV-2-9-B	100pcs/box
			"-" Pre-slit	NKCSV-2-9-B-P	100pcs/box
		Blue PTFE/White Silicone	non-slit	NKCSV-2-9-F	100pcs/box
			"-" Pre-slit	NKCSV-2-9-F-P	100pcs/box
ND9 Sampler Vials Kit (assembled) Amber vial +blue cap	2ml Amber vials w/writing spot	White PTFE/Red Silicone	non-slit	NKASV-2-9-A	100pcs/box
			"-" Pre-slit	NKASV-2-9-A-P	100pcs/box
		Red PTFE/White Silicone	non-slit	NKASV-2-9-B	100pcs/box
			"-" Pre-slit	NKASV-2-9-B-P	100pcs/box
		Blue PTFE/White Silicone	non-slit	NKASV-2-9-F	100pcs/box
			"-" Pre-slit	NKASV-2-9-F-P	100pcs/box

LC/MS and GC/MS Certified Vial Kits

Certification Support: Every Piece Meets Professional Standards

Full batches have passed USP Class VI Biocompatibility Certification, ISO 9001 Quality Management System Certification and RoHS Environmental Certification. Complete batch test reports (including extractable testing and particle content analysis) are provided.

Among them, USP Class VI Certification ensures no harmful substance migration when in contact with food and pharmaceuticals, fully meeting the sample processing and storage requirements under GLP standards.



Description	Septum	Slit Type	Cat. No.	Pack
Certified Kits, Pre-sealed, 2mL Clear Screw Vial w/Marking Spot+Blue Cap	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKCASV-2-9-A	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKCASV-2-9-A-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKCASV-2-9-B	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKCASV-2-9-B-P	100pcs/bag
Certified Kits, Pre-sealed, 2mL Amber Screw Vial w/Marking Spot+Blue Cap	White PTFE/Red Silicone, 9*1.0mm	non-slit	NKCCSV-2-9-A	100pcs/bag
	White PTFE/Red Silicone, 9*1.0mm	"-" Pre-slit	NKCCSV-2-9-A-P	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	non-slit	NKCCSV-2-9-B	100pcs/bag
	Red PTFE/White Silicone, 9*1.0mm	"-" Pre-slit	NKCCSV-2-9-B-P	100pcs/bag

PP Screw Closures ND10 with Septa

Thread Standard: 10mm outer diameter of the bottle mouth paired with Type 425 thread; adopts a screw-thread sealing design for reliable sealing and reusability.

Easy Identification: Equipped with a writing area and clear volume markings at 0.5, 1.0 and 1.5mL, facilitating sample labeling and quantitative sampling.

Excellent Compatibility: Boasts high dimensional accuracy and is compatible with autosamplers of mainstream instruments including Agilent, Shimadzu, Waters for HPLC, GC, LC-MS applications.



NKSSC-10-A NKSSC-10-A-P NKSSC-10-B NKSSC-10-B-P NKSV-2-10 NKSV-2-10-W NKSV-2-10-A NKSV-2-10-A-W
NKSV-2-10-U NKSV-2-10-W-U NKSV-2-10-A-U NKSV-2-10-A-W-U

Description	Septum	Slit Type	Cat. No.	Pack
10mm Black Screw Polypropylene Cap, 7mm Centre Hole	White PTFE/Red Silicone Septa	non-slit	NKSSC-10-F	100pcs/bag
	White PTFE/Red Silicone Septa	"-" Pre-slit	NKSSC-10-S	100pcs/bag
	Red PTFE/White Silicone Septa	non-slit	NKSSC-10-A	100pcs/bag
	Red PTFE/White Silicone Septa	"-" Pre-slit	NKSSC-10-B	100pcs/bag

Screw Neck Vials ND10

Description	Material	Writing area	Cat. No.	Pack
2ml, 10-425 Screw Clear Vial, 11.6*32mm	Boro 7.0	NO	NKSV-2-10	100pcs/box
	Boro 5.0	NO	NKSV-2-10-U	100pcs/box
	Boro 7.0	YES	NKSV-2-10-W	100pcs/box
	Boro 5.0	YES	NKSV-2-10-W-U	100pcs/box
2ml, 10-425 Screw Amber Vial, 11.6*32mm	Boro 7.0	NO	NKSV-2-10	100pcs/box
	Boro 5.0	NO	NKSV-2-10-A-U	100pcs/box
	Boro 7.0	YES	NKSV-2-10-A-W	100pcs/box
	Boro 5.0	YES	NKSV-2-10-A-W-U	100pcs/box

Silver Crimp-top Aluminum Cap

The 11mm crimp vial septa is the mainstream choice for micro-sample analysis in laboratories. Precisely compatible with 1mL/2mL borosilicate 3.3 glass vials, it has become a standardized accessory for routine GC/LC testing and trace analysis, thanks to its dual advantages of rigid sealing and precise puncture.



NKSSCC-11-A



NKSSCC-11-B



NKSSCC-11-C



NKSSCC-M-11-A



NKSSCC-M-11-B



NKSSCC-M-11-R

Description	Septum	Cat. No.	Pack
ND11 Silver Crimp-top Aluminum Cap, with central hole 5.5mm	White PTFE/Red Silicone, 11*1mm	NKSSCC-11-A	100pcs/bag
	Red PTFE/White Silicone, 11*1mm	NKSSCC-11-B	100pcs/bag
	Red PTFE/Sil/Red PTFE	NKSSCC-11-C	100pcs/bag
ND11 Silver Crimp-top Aluminum Cap, 5.5mm Centre Hole (new)	White PTFE/Red Silicone, 11*1mm	NKSSCC-M-11-A	100pcs/bag
	Red PTFE/White Silicone, 11*1mm	NKSSCC-M-11-B	100pcs/bag
	Transparent PTFE/Orange Natural Rubber Septa, 11*1mm	NKSSCC-M-11-R	100pcs/bag

Crimp Neck Vials ND11

Description	Material	Writing area	Cat. No.	Pack
2ML ND11 Clear Glass Crimp Neck Vial, 11.6*32mm	Boro 7.0	NO	NKCV-2-11	100pcs/box
	Boro 7.0	NO	NKCV-2-11-W	100pcs/box
	Boro 5.0	YES	NKCV-2-11-U	100pcs/box
	Boro 5.0	YES	NKCV-2-11-W-U	100pcs/box
2ML ND11 Amber Glass Crimp Top Vial, 11.6*32mm	Boro 7.0	NO	NKCV-2-11-A	100pcs/box
	Boro 7.0	NO	NKCV-2-11-A-W	100pcs/box
	Boro 5.0	YES	NKCV-2-11-A-U	100pcs/box
	Boro 5.0	YES	NKCV-2-11-A-W-U	100pcs/box

Micro-Vials ND11 for Small Sample Volumes (High Recovery Vials)



NKCV-03-11



NKCV-11-11



NKCV-13-11



NKCV-15-11



NKCV-03-11-A



NKCV-11-11-A



NKCV-13-11-A



NKCV-15-11-A

Description	Volume	Cat. No.	Pack
ND11 Clear Crimp Glass Vial High recovery, 11.6x 32mm	0.3ml Integrated w/Micro-Insert Base Bonded	NKCV-03-11	100pcs/box
	1.1 mL vial, 25 μ L Reservoir	NKCV-11-11	100pcs/box
	1.3 mL, 30 μ L Reservoir	NKCV-13-11	100pcs/box
	1.5ml Champagne Shape base	NKCV-15-11	100pcs/box
ND11 Amber Crimp Glass Vial High recovery, 11.6x 32mm	0.3ml Integrated w/Micro-Insert Base Bonded	NKCV-03-11-A	100pcs/box
	1.1 mL vial, 25 μ L Reservoir	NKCV-11-11-A	100pcs/box
	1.3 mL, 30 μ L Reservoir	NKCV-13-11-A	100pcs/box
	1.5ml Champagne Shape base	NKCV-15-11-A	100pcs/box

PP Snap Closures ND11 with Septa

Snap top vials and snap caps are a sample containment solution for nonvolatile liquid samples.

Snap top vials enable the user to hold both the pipette and snap cap in one hand and the snap cap vial in the other when pipetting, making them more convenient than other vial types.



NKSSNC-11-A



NKSSNC-11-B



NKSSNC-11-F



NKSSNC-11-A-P





NKSSNC-11-B-P



NKSSNC-11-F-P

Description	Septum	Slit Type	Cat. No.	Pack
11mm PE Snap Cap, 6mm Centre Hole	White PTFE/Red Silicone, 11*1mm	non-slit	NKSSNC-11-A	100pcs/bag
	Red PTFE/White Silicone, 11*1mm	non-slit	NKSSNC-11-B	100pcs/bag
	Blue PTFE/White Silicone, 11*1mm	non-slit	NKSSNC-11-F	100pcs/bag
	White PTFE/Red Silicone, 11*1mm	"-" Pre-slit	NKSSNC-11-A-P	100pcs/bag
	Red PTFE/White Silicone, 11*1mm	"-" Pre-slit	NKSSNC-11-B-P	100pcs/bag
	Blue PTFE/White Silicone, 11*1mm	"-" Pre-slit	NKSSNC-11-F-P	100pcs/bag

Snap Ring Vials ND11

Description		Material	Writing area	Cat. No.	Pack
2mL Clear Snap Top Glass Vial, 11.6*32mm		Boro 7.0	NO	NKNV-2-11	100pcs/box
		Boro 7.0	YES	NKNV-2-11-W	100pcs/box
		Boro 5.0	NO	NKNV-2-11-U	100pcs/box
		Boro 5.0	YES	NKNV-2-11-W-U	100pcs/box
2mL Amber Snap Top Glass Vial, 11.6*32mm		Boro 7.0	NO	NKNV-2-11-A	100pcs/box
		Boro 7.0	YES	NKNV-2-11-A-W	100pcs/box
		Boro 5.0	NO	NKNV-2-11-A-U	100pcs/box
		Boro 5.0	YES	NKNV-2-11-A-W-U	100pcs/box

PP Screw Closures ND13 with Septa

Compliant with industry size standards, compatible with automatic samplers from brands such as Agilent, Waters, and Shimadzu. Ready-to-use, enhancing experimental workflow efficiency.



NKSSC-13-A



NKSSC-13-B



NKSSC-13-A-C



NKSSC-13-B-C



NKSSC-13-H-C



NKSC-13-S

Description	Septum	Cat. No.	Pack
13mm Black Screw Polypropylene Cap with Septa, 8.5mm Centre Hole	White PTFE/Red Silicone, 13*1.5mm	NKSSC-13-A	100pcs/bag
	Red PTFE/White Silicone, 13*1.5mm	NKSSC-13-B	100pcs/bag
13mm Screw Polypropylene Black Cap, close top	White PTFE/Red Silicone, 13*1.5mm	NKSSC-13-A-C	100pcs/bag
	Red PTFE/White Silicone, 13*1.5mm	NKSSC-13-B-C	100pcs/bag
	Nature PTFE/Nature Silicone, 1.5mm thickness	NKSSC-13-H-C	100pcs/bag
13-425 Black Screw cap with conical plug		NKSC-13-S	100pcs/bag

Description		Material	Writing area	Cat. No.	Pack
4mL Clear Screw Neck Vial, 13-425 Thread, 14.5*45mm		Boro 7.0	NO	NKSV-4-13	100pcs/box
		Boro 7.0	YES	NKSV-4-13-W	100pcs/box
		Boro 5.0	NO	NKSV-4-13-U	100pcs/box
		Boro 5.0	YES	NKSV-4-13-W-U	100pcs/box
4mL Amber Screw Neck Vial, 13-425 Thread, 14.5*45mm		Boro 7.0	NO	NKSV-4-13-A	100pcs/box
		Boro 7.0	YES	NKSV-4-13-A-W	100pcs/box
		Boro 5.0	NO	NKSV-4-13-A-U	100pcs/box
		Boro 5.0	YES	NKSV-4-13-A-W-U	100pcs/box

Screw Neck ND15

PP Screw Closures ND15 with Septa




Description	Septum	Cat. No.	Pack
15-425 Black Screw Polypropylene Cap, 9mm Centre Hole	White PTFE/Red Silicone, 1.5mm thickness	NKSSC-15-A	100pcs/bag
	Red PTFE/White Silicone, 1.5mm thickness	NKSSC-15-B	100pcs/bag
15-425 Black Screw Polypropylene Cap, closed-top	White PTFE/Red Silicone, 1.5mm thickness	NKSSC-15-A-C	100pcs/bag
	Red PTFE/White Silicone, 1.5mm thickness	NKSSC-15-B-C	100pcs/bag

Screw Neck Vials ND15

Description	Details	Cat. No.	Pack
8ML ND15 Glass Screw Vial, 15-425 Thread	Clear, 16.6*60mm	NKSV-8-15	100pcs/box
	Amber, 16.6*60mm	NKSV-8-15-A	100pcs/box
12ML ND15 Glass Screw Vial, 15-425 Thread	Clear, 16.6*60mm	NKSV-12-15	100pcs/box
	Amber, 16.6*60mm	NKSV-12-15-A	100pcs/box

Magnetic Screw Closures ND18 with Septa

Precision thread screw seals ND18 are metal caps with fitted septa, compatible with SPME and headspace applications. The design ensures sufficient surface for magnetic transport while maintaining a flat, secure fit. Easier to dispose than crimp caps. These closed tops caps are ideal for sample storage when paired with ND18 vials.



Description	Details	Cat. No.	Pack
18mm Silver Magnetic Precision Screw Metal Cap, 8mm Centre Hole, with Septa	Blue PTFE/White Silicone, 17.5*1.5mm	NKSSC-18M-F	100pcs/bag
	White PTFE/Light Blue Silicone, 17.5*1.3mm	NKSSC-18M-L	100pcs/bag
	White PTFE/Transparent Blue Silicone, 17.5*1.5mm, ultra clean	NKSSC-18M-T	100pcs/bag
	PFAS suitable, Transparent Silicone / Transparent Polypropylene, 17.5X1.5mm	NKSSC-18M-TT	100pcs/bag
18mm Magnetic Precision Screw Metal Cap, close top	Blue PTFE/White Silicone Septa 17.5X1.5mm	NKSSC-18M-F-C	100pcs/bag

Precision Thread Screw Neck Vials ND18



Description	Color	Material	Cat. No.	Pack
10ml Round Bottom, 22.5*46mm	Clear	Boro 7.0	NKSV-10R	100pcs/box
	Clear	Boro 5.0	NKSV-10R-U	100pcs/box
10ml Round Bottom, 22.5*46mm	Amber	Boro 7.0	NKSV-10R-A	100pcs/box
	Amber	Boro 5.0	NKSV-10R-A-U	100pcs/box
20ml Round Bottom, 22.5*75.5mm	Clear	Boro 7.0	NKSV-20R	100pcs/box
	Clear	Boro 5.0	NKSV-20R-U	100pcs/box
20ml Round Bottom, 22.5*75.5mm	Amber	Boro 7.0	NKSV-20R-A	100pcs/box
	Amber	Boro 5.0	NKSV-20R-A-U	100pcs/box

Aluminum Crimp Closures and Magnetic Crimp Closures ND20 with Septa

Dual Closure Options for Diverse Scenarios


Aluminum Crimp Closures: Corrosion-resistant aluminum alloy, pressure-resistant and airtight, prevents sample volatilization and leakage, suitable for manual/semiautomatic workflows.
 Magnetic Crimp Closures: Built-in magnetic components, compatible with automatic samplers from CTC, Metrohm, Shimadzu, Agilent, etc., enhancing high-throughput efficiency.

ND20 Standard Size, Universal Compatibility

Compliant with industry standards, 20mm diameter fits ND20 series vials, compatible with HPLC/GC/MS and other mainstream instruments, ready-to-use.

Composite Septa, Enhanced Sealing

PTFE/silicone composite material, temperature resistance: -60 C ~200 C, corrosion-resistant and non-contaminating, leak-proof sealing, protects autosampler needles.



NKCC-20

NKSSCC-20-G

NKSSCC-20-H

NKSSCC-20-L

NKSSCC-20-T-U

NKSSCC-20-W

NKSSCC-20M-E

NKSSCC-20M-S

NKSSCC-20M-T-U

NKSSCC-20M-R

Headspace Crimp Neck ND20

Description	Septum	Cat. No.	Pack
Crimp Cap ND20, Empty	No septa	NKCC-20	100pcs/box
20mm Crimp-top Aluminum Cap w/septa 9.5mm Centre Hole	Beige PTFE/White Silicon, 20*3.0mm	NKSSCC-20-G	100pcs/box
	Nature PTFE/Nature Silicone, 20*3.0mm	NKSSCC-20-H	100pcs/box
	White PTFE/Light Blue Silicone, 20*3mm	NKSSCC-20-L	100pcs/box
	Transparent PTFE/Transparent Blue Silicone, 20*3.0mm, Ultra Clean	NKSSCC-20-T-U	100pcs/box
	White PTFE/White Silicone, 20*3mm	NKSSCC-20-W	100pcs/box

Description	Septum	Cat. No.	Pack
20mm Crimp-top Blue Magnetic Aluminum Cap, 8mm Centre Hole	White PTFE/Blue Silicone Septa 20X3.0mm	NKSSCC-20M-E	100pcs/box
	Nature PTFE/White Silicone Septa 20X3.0mm	NKSSCC-20M-S	100pcs/box
	Transparent PTFE/Blue transparent Silicone Septa 20X3.0mm, ultra clean	NKSSCC-20M-T-U	100pcs/box
20mm Crimp-top Red Magnetic Aluminum Cap, 8mm Centre Hole	Nature PTFE/White Silicone Septa 20X3.0mm	NKSSCC-20M-R	100pcs/box



NKSSCC-20S-G



NKSSCC-20S-H



NKSSCC-20S-R



NKSSCC-20S-W

Description	Septum	Cat. No.	Pack
20mm Pressure Release Aluminum Crimp Cap With Septa, 8mm Centre Hole	Beige PTFE/White Silicon, 20*3.0mm	NKSSCC-20S-G	100pcs/box
	Nature PTFE/Nature Silicone, 20*3.0mm	NKSSCC-20S-H	100pcs/box
	Natural PTFE/White Silicone, 20*3mm	NKSSCC-20S-R	100pcs/box
	White PTFE/White Silicone, 20*3mm	NKSSCC-20S-W	100pcs/box



NKCV-6F
NKCV-6F-U



NKCV-6F-1
NKCV-6F-1-U



NKCV-10F
NKCV-10F-U



NKCV-10F-A
NKCV-10F-A-U



NKCV-10R
NKCV-10R-U



NKCV-10R-A
NKCV-10R-A-U

Description	Material	Shape	Cat. No.	Pack
6ml/7ml Clear Crimp Headspace Vial, 22.5*38mm	Boro 7.0	Flat Bottom	NKCV-6F	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-6-F-U	100pcs/box
6mL Clear Crimp Headspace Vial, 21.7*38mm	Boro 7.0	Flat Bottom	NKCV-6-F-1	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-6-F-1-U	100pcs/box

Headspace Crimp Neck ND20

Description	Material	Shape	Cat. No.	Pack
ND20, 10 mL Crimp Headspace Clear Vial 22.5*46mm	Boro 7.0	Flat Bottom	NKCV-10F	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-10F-U	100pcs/box
ND20, 10 mL Crimp Headspace Amber Vial 22.5*46mm	Boro 7.0	Flat Bottom	NKCV-10F-A	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-10F-A-U	100pcs/box
ND20, 10 mL Crimp Headspace Clear Vial 22.5*46mm	Boro 7.0	Round bottom	NKCV-10R	100pcs/box
	Boro 5.0	Round bottom	NKCV-10R-U	100pcs/box
ND20, 10 mL Crimp Headspace Amber Vial 22.5*46mm	Boro 7.0	Round bottom	NKCV-10R-A	100pcs/box
	Boro 5.0	Round bottom	NKCV-10R-A-U	100pcs/box



NKCV-20F NKCV-20F-75 NKCV-20F-A NKCV-20F-75-A NKCV-20R NKCV-20R-75 NKCV-20R-A NKCV-20R-75-A
 NKCV-20F-U NKCV-20F-75-U NKCV-20F-A-U NKCV-20F-75-A-U NKCV-20R-U NKCV-20R-75-U NKCV-20R-A-U NKCV-20R-75-A-U

Description	Material	Shape	Cat. No.	Pack
ND20, 20 mL Crimp Headspace Clear Vial 22.5*75mm	Boro 7.0	Flat Bottom	NKCV-20F	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-20F-U	100pcs/box
ND20, 20mL Crimp Headspace Clear Vial 22.75*75mm	Boro 7.0	Flat Bottom	NKCV-20F-75	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-20F-75-U	100pcs/box
ND20, 20 mL Crimp Headspace Amber Vial 22.5*75mm	Boro 7.0	Flat Bottom	NKCV-20F-A	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-20F-A-U	100pcs/box
ND20, 20 mL Crimp Headspace Amber Vial 22.75*75mm	Boro 7.0	Flat Bottom	NKCV-20F-75-A	100pcs/box
	Boro 5.0	Flat Bottom	NKCV-20F-75-A-U	100pcs/box
ND20, 20 mL Crimp Headspace Clear Vial 22.5*75mm	Boro 7.0	Round bottom	NKCV-20R	100pcs/box
	Boro 5.0	Round bottom	NKCV-20R-U	100pcs/box
ND20, 20 mL Crimp Headspace Clear Vial 22.75*75mm	Boro 7.0	Round bottom	NKCV-20R-75	100pcs/box
	Boro 5.0	Round bottom	NKCV-20R-75-U	100pcs/box
ND20, 20 mL Crimp Headspace Amber Vial 22.5*75mm	Boro 7.0	Round bottom	NKCV-20R-A	100pcs/box
	Boro 5.0	Round bottom	NKCV-20R-A-U	100pcs/box
ND20, 20mL Crimp Headspace Amber Vial 22.75*75mm	Boro 7.0	Round bottom	NKCV-20R-75-A	100pcs/box
	Boro 5.0	Round bottom	NKCV-20R-75-A-U	100pcs/box

PP Screw Closures ND24



Description	Septum	Cat. No.	Pack
24-400 White Screw PP cap with 15mm central hole		NKSC-22	100pcs//bag
24-400 Black PP Screw cap with conical plug		NKSC-24-S	100pcs//bag
24-400 White Screw PP Cap, 15mm Centre Hole with Septae	Beige PTFE/White Silicon, 22*3.0mm	NKSSC-22-G	100pcs//bag
	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSC-22-H	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22-S	100pcs//bag
24-400 Black Screw PP Cap with 15mm Centre Hole	Beige PTFE/White Silicon, 22*3.0mm	NKSSC-22-G-1	100pcs//bag
	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSC-22-H-1	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22-S-1	100pcs//bag



NKSSC-22-X-C NKSSC-22-S1-C NKSSC-22-H1-C NKSSC-22-G-1-C NKSSC-22-H-1-C NKSSC-22-S-1-C NKSSC-22-X-1-C

Description	Septum	Cat. No.	Pack
24-400 White Screw PP Cap, Closed-top	Beige PTFE/White Silicon, 22*3.0mm	NKSSC-22-G-C	100pcs//bag
	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSC-22-H-C	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22-S-C	100pcs//bag
24-400 White Screw PP Cap, Closed-top	Natural PTFE/Red silicone, 22*3mm	NKSSC-22-X-C	100pcs//bag
	Natural PTFE/White Silicone, 22*1.5mm	NKSSC-22-S1-C	100pcs//bag
	Nature PTFE/Nature Silicone, 1.5mm thickness	NKSSC-22-H1-C	100pcs//bag
24-400 Black Screw PP Cap, Closed-top	Beige PTFE/White Silicon, 22*3.0mm	NKSSC-22-G-1-C	100pcs//bag
	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSC-22-H-1-C	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22-S-1-C	100pcs//bag
	Natural PTFE/Red Silicon, 22*3mm	NKSSC-22-X-1-C	100pcs//bag

UltraBond Closures ND24



Description	Septum	Cat. No.	Pack
Bonded, 24-400 White Screw PP Cap with Septa, 15mm Centre Hole	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSC-22B-H	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSC-22B-S	100pcs//bag
Bonded, 24-400 White Screw PP Cap with Septa, 16mm Centre Hole	Nature PTFE/Nature Silicone, 22*3.0mm	NKSSBC-22B-H	100pcs//bag
	Nature PTFE/White Silicone, 22*3.0mm	NKSSBC-22B-S	100pcs//bag
	Natural PTFE/White Silicone, 22*1.5mm	NKSSBC-22B-S-15	100pcs//bag

Screw Neck Vials ND24 for TOC / VOC Analysis and Sample Storage

NKSV-20-E
NKSV-20-E-UNKSV-20-E-A
NKSV-20-E-A-UNKSV-30-E
NKSV-30-E-UNKSV-30-E-A
NKSV-30-E-A-UNKSV-40-E
NKSV-40-E-UNKSV-40-E-A
NKSV-40-E-A-UNKSV-60-E
NKSV-60-E-UNKSV-60-E-A
NKSV-60-E-A-U

Description	Material	Color	Cat. No.	Pack
ND24 20ml Screw Neck Vial, 27.5x57mm	Boro 7.0	Clear	NKSV-20-E	100pcs/box
	Boro 5.0	Clear	NKSV-20-E-U	100pcs/box
	Boro 7.0	Amber	NKSV-20-E-A	100pcs/box
	Boro 5.0	Amber	NKSV-20-E-A-U	100pcs/box
ND24 30ml Screw Neck Vial, 27.5x75mm	Boro 7.0	Clear	NKSV-30-E	100pcs/box
	Boro 5.0	Amber	NKSV-30-E-U	100pcs/box
	Boro 7.0	Clear	NKSV-30-E-A	100pcs/box
	Boro 5.0	Amber	NKSV-30-E-A-U	100pcs/box
ND24 40ml Screw Neck Vial, 27.5x95mm	Boro 7.0	Clear	NKSV-40-E	100pcs/box
	Boro 5.0	Amber	NKSV-40-E-U	100pcs/box
	Boro 7.0	Clear	NKSV-40-E-A	100pcs/box
	Boro 5.0	Amber	NKSV-40-E-A-U	100pcs/box
ND24 60ml Screw Neck Vial, 27.5x140mm	Boro 7.0	Clear	NKSV-60-E	100pcs/box
	Boro 5.0	Amber	NKSV-60-E-U	100pcs/box
	Boro 7.0	Clear	NKSV-60-E-A	100pcs/box
	Boro 5.0	Amber	NKSV-60-E-A-U	100pcs/box

TOC Sets: Pre-Cleaned and Pre-Sealed Vials with Closures ND24

A dedicated laboratory consumable engineered for precise determination of Total Organic Carbon (TOC) and Volatile Organic Compounds (VOCs), serving as a reliable guarantee throughout the entire process of sample collection, storage, and transportation.

Premium Materials to Eliminate Contamination

The bottle caps and septa are made of PTFE/silicone composite material, featuring excellent corrosion resistance and the ability to withstand repeated punctures by auto-sampler needles. With superior sealing performance, they ensure zero sample leakage.

Standard Specifications with Wide Compatibility

Available in multiple capacities including 20mL, 30mL, and 40mL (40mL as the industry's mainstream specification) to meet the needs of various testing scenarios. The bottle mouth adopts a standard 24-400 thread design with an outer diameter of approximately 27.5mm, enabling seamless compatibility with a variety of TOC analyzers and purge and trap systems, thus significantly reducing equipment adaptation costs.



Description	Dust Proof Cap	Central Hole	Volume	Color	Cat. No.	Pack
ND24 Screw Vial + White PP Cap with 3mm Nature PTFE/ White Sil Septa	NO	15mm	20ml	Clear	NKTOC-20C	100pcs/box
	NO	15mm	20ml	Amber	NKTOC-20A	100pcs/box
	YES	15mm	20ml	Clear	NKTOC-20C-DC	100pcs/box
	YES	15mm	20ml	Amber	NKTOC-20A-DC	100pcs/box
	NO	15mm	40ml	Clear	NKTOC-40C	100pcs/box
	NO	15mm	40ml	Amber	NKTOC-40A	100pcs/box
	YES	15mm	40ml	Clear	NKTOC-40C-DC	100pcs/box
	YES	15mm	40ml	Amber	NKTOC-40A-DC	100pcs/box
	NO	15mm	60ml	Clear	NKTOC-60C	100pcs/box
	NO	15mm	60ml	Amber	NKTOC-60A	100pcs/box
	YES	15mm	60ml	Clear	NKTOC-60C-DC	100pcs/box
	YES	15mm	60ml	Amber	NKTOC-60A-DC	100pcs/box
ND24 Screw Vial + White PP Cap with 1.5mm Nature PTFE/ White Sil Septa	YES	15mm	40ml	Clear	NKTOC-40C-DC-15	100pcs/box
	YES	15mm	40ml	Amber	NKTOC-40A-DC-15	100pcs/box
ND24 Screw Vial + White PP bonded Cap with 3mm Nature PTFE/White Sil Septa	YES	15mm	40ml	Clear	NKTOC-40CB-DC	100pcs/box
	YES	15mm	40ml	Amber	NKTOC-40AB-DC	100pcs/box
ND24 Screw Vial + White PP bonded Cap with 1.5mm Nature PTFE/White Sil Septa	YES	15mm	40ml	Clear	NKTOC-40CB-DC	100pcs/box
	YES	15mm	40ml	Amber	NKTOC-40AB-DC	100pcs/box

COD Vials and Closures



NKPCC-N-PR-1



NKPCC-N-PR-1A



NKPCC-N-PR-2



NKPCC-N-PR-2A



NKPPC-CD-W-PI-1



NKPPC-W-PR-4

Description	Septum	Temperature	Cat. No.	Pack
COD Nature PC 15-415 Cap, Top Feed	Nature PTFE/RU/Nature PTFE	≤175°C	NKPCC-N-PR-1	100pcs//bag
	Nature PTFE/RU/Nature PTFE, With 120° Logo	≤175°C	NKPCC-N-PR-1A	100pcs//bag
	Natural PTFE/RU/Blue PTFE	≤175°C	NKPCC-N-PR-2	100pcs//bag
	Natural PTFE/RU/Blue PTFE, With 120° Logo	≤175°C	NKPCC-N-PR-2A	100pcs//bag
COD White PP 15-415 Cap	Nature PTFE/SIL/Nature PTFE Low bleed Septa	≤150°C	NKPPC-CD-W-PI-1	100pcs//bag
COD White PP GL14 Cap	With Gray PTFE/RU/Gray PTFE	≤150°C	NKPPC-W-PR-4	100pcs//bag



NKTT-COD-30F-U-100



NKTT-COD-30R-U-100



NKTT-COD-30R-U-125



NKTT-COD-30F-U-90



NKTT-COD-30F-U-100-N

Description	Shape	Cat. No.	Pack
COD Glass Tube Clear, 15-415, 10mL, 16x100mm, HC type	Flat Bottom	NKTT-COD-30F-U-100	100pcs/box
	Round bottom	NKTT-COD-30R-U-100	100pcs/box
COD Glass Tube Clear, 15-415, 10mL, 16x125mm, HC type	Round bottom	NKTT-COD-30R-U-125	100pcs/box
COD Glass Tube Clear, GL14, 9mL, 16x90mm, LV type	Flat Bottom	NKTT-COD-30F-U-90	100pcs/box
COD Glass Tube Clear, GL14, 10mL, 16X97mm, HN type	Flat Bottom	NKTT-COD-30F-U-100-N	100pcs/box

Shell Vials and Closures





Designed for high-throughput lab routine analysis, Shell Vials & Closures feature a neckless straight-wall design and push-in sealing for fast operation. Cost-effective and easy to use, they're ideal for routine HPLC, sample pretreatment, non-volatile sample storage, and compatible with mainstream auto-injectors.

NKSU-1

1ml 8x40 Flat Base Clear Glass Shell Vial with PE Plug

Sample Cell Vials with Closures

Made of high borosilicate glass, ensuring accurate and repeatable turbidity measurements for water quality and environmental monitoring.

Description	Septum	Marking	Cat. No.	Pack
Sample Cell tube 20ml, 25x60mm, assembled with PF cap 	PE foam/PTFE septa	no marking	NKSCT-2560	100pcs/box
	PE foam/PTFE septa	with 10ml marking	NKSCT-2560-10G	100pcs/box
	PE foam/PTFE septa	with 15ml marking	NKSCT-2560-15G	100pcs/box
Sample Cell tube 30ml, 25x95mm, assembled with PF cap 	pure silicone septa	no marking	NKSCT-2560-S	100pcs/box
	pure silicone septa	with 10ml marking	NKSCT-2560-10G-S	100pcs/box
	pure silicone septa	with 15ml marking	NKSCT-2560-15G-S	100pcs/box
	PE foam/PTFE septa	no marking	NKSCT-2595	100pcs/box
	pure silicone septa	no marking	NKSCT-2595-S	100pcs/box

Scintillation Vials with Closures

- For scintillation counting and standard storage applications.
- Packaged in convenient utility trays, for easy sample storage.
- Use vials with PP cap and metal foil / pulp liner for applications involving hydrocarbons, oils, ketones and alcohol.
- Urea caps with PE cone liner offer superior torque retention and excellent sealing characteristics.

NKSCT-20-VC	20mL Clear Glass Scintillation Vial, 27.5x57.5mm, 22-400 Screw Thread with PP cap and PE septa & Aluminum foil gasket
NKSCT-20-VC-1	20mL Clear Glass Scintillation Vial, 27.5x57.5mm, 22-400 Screw Thread with Urea Cap & Cone PE Plug



NMR Tubes

- Outer diameter: 4.93±0.03mm; Length: 7inch (178mm) – perfect rotor compatibility
- Wall thickness: 0.43mm (maximizes strength, minimizes cracking)
- Material: Borosilicate 3.3 glass (ASTM-compliant, acid/alkali-resistant)
- Centrifuge runout <0.08mm (maintains strict cylindrical shape during rotation)
- Average convexity: 60µm (optimizes resolution for <250 molecular weight samples at up to 600MHz)
- Compatible rotors: Wilmad STB-5, Brooke Z42516, JEOL 800261275



NKNMR-5-B	NMR tube, length 178mm (7inch), dia. 4.93mm, boro material, with colorful caps	100 pcs/box
NKNMR-C	Colorful caps for 4.93mm I.D. dia NMR tube, normal Type	100 pcs/box
NKNMR-C1	Colorful caps for 4.93mm I.D. dia NMR tube, Thread Type	100 pcs/box

Blue screw cap glass reagent bottle


Reagent bottle made by auto-linear machine from germany, Which protect the weight, temperature, pressure be uinfrom. And the bottle has no leak, thick fineness. It consists of special component that has stable chemistry, and not go moule.



Clear scale



Strong sealing performance



Thickened material





Description	Volume	Cat. No.	Pack
Glass blue screw cap reagent bottle, clear, soda lime glass	25ml	1407-25	288pcs/ctn
	50ml	1407-50	144pcs/ctn
	100ml	1407-100	96pcs/ctn
	250ml	1407-250	70pcs/ctn
	500ml	1407-500	24pcs/ctn
	1000ml	1407-1000	15pcs/ctn
	2000ml	1407-2000	12pcs/ctn
	5000ml	1407-5000	6pcs/ctn
Glass blue screw cap reagent bottle, amber, soda lime glass	100ml	1407/1-100	96pcs/ctn
	250ml	1407/1-250	70pcs/ctn
	500ml	1407/1-500	24pcs/ctn
	1000ml	1407/1-1000	15pcs/ctn

Description	Volume	Cat. No.	Pack
Glass blue screw cap reagent bottle, Clear, boro 3.3	50ml	1408-50	162pcs/ctn
	100ml	1408-100	90pcs/ctn
	250ml	1408-250	56pcs/ctn
	500ml	1408-500	35pcs/ctn
	1000ml	1408-1000	24pcs/ctn
Glass blue screw cap reagent bottle, Clear, boro3.3, GL80	1000ml	1409-1000	24pcs/ctn
	10000ml	1409-10000	4pcs/ctn
	2000ml	1409-2000	12pcs/ctn
	5000ml	1409-5000	6pcs/ctn
	20000ml	1409-20000	1pcs/ctn

HPLC Safety Caps & Waste Container Caps

Solvent Safety Cap is equipped with an original one-way air intake valve design. Each solvent tube is fitted with a dedicated pipeline connector, ensuring secure fixation and tight sealing of the pipeline. When the instrument is in operation, the pump extracts solvent to generate negative pressure, which triggers the one-way valve to open automatically, allowing filtered clean air to enter the solvent bottle. After the instrument shuts down, the one-way valve closes immediately, effectively blocking solvent volatilization and moisture absorption, and fully safeguarding experimental safety and solvent purity throughout the process.





GL45 solvent Safety Cap



NKSSC0145	NKSSC0245	NKSSC0345	NKSSC0445	NKSSC0545
GL45 Single port solvent Safety Cap Including: 1x Air Valve+1x 3.2mm pipe joint	GL45 2 ports solvent Safety Cap Including: 1x Air Valve+2x 3.2mm pipe joint+1x plug	GL45 3 ports solvent Safety Cap Including: 1x Air Valve+3x 3.2mm pipe joint+2x plugs	GL45 4 ports solvent Safety Cap Including: 1x Air Valve+4x 3.2mm pipe joint+3x plugs	GL45 5 ports solvent Safety Cap Including: 1x Air Valve+5x 3.2mm pipe joint+4x plugs

Waste Liquid Safety Cap,

Height: 55mm, O.D.: 72mm, S60 Mouth, Joint Thread 9.2mm




NKWLSC0160	NKWLSC0160D	NKWLSC0360	NKWLSC0460	NKWLSC0560	NKWLSC0660
Single port, with 1x soft waste pipe joint+1x Y joint	Single port, with Alarm D, with 1x soft waste pipe joint+1x Y joint	3 ports, 1x soft waste pipe joint+2x 3.2mm pipe joint+ 2x plug+1x Y joint	4 ports, with Alarm D, with 1x soft waste pipe joint+3x 3.2mm pipe joint+3x plug+ 1x Y joint	5 ports, 1x soft waste pipe joint+4x 3.2mm pipe joint+ 4x plug+1x Y joint	6 ports, 1x soft waste pipe joint+5x 3.2mm pipe joint+ 5x plug+1x Y joint

Air Valve (1 X 1/4-28 Fitting)				Waste Liquid Bucket, White, HDPE, S60 mouth			
NKAV00A	NKAV00B	NKAV00C	NKEF014	NKWLB0160-10	NKWLB0160-20	NKWLB0160-20S	NKWLB0160-30
Air Valve (1 X 1/4-28 Fitting) A Style	Air Valve (1 X 1/4-28 Fitting) B Style with Writing Label	Air Valve (1 X 1/4-28 Fitting) C Style with Timing Label	GL14-65g Exhaust Filter, 13.5cm length, 6cm diameter	10L, 0.55kg, LxWxH : 23*20*28cm	20L, 1.3kg, LxWxH : 27*27*39cm	20L short, 1.3kg LxWxH: 32*29*30cm	30L, 1.7kg, LxWxH : 29*32*43cm

Racks for Vials

Designed specifically for laboratory vials (including HPLC/GC dedicated vials), it combines stability, compatibility and convenience to facilitate standardized storage, safe protection and efficient access of samples. It is an essential consumable for chromatographic analysis, sample pre-treatment, routine experiments and low-temperature storage scenarios.

Featuring standard-preset hole diameters and slot designs, it precisely fits common vial sizes such as 1.5mL, 2mL, 5mL and 10mL, and is compatible with screw-top vials, snap-top vials and chromatography-specific vials (e.g., HPLC injection vials).



Name	Description	Cat. No.	Pack
Vials Rack	6-20ml Vials Rack 40 Wells, Well Dia 24mm, for Crimp Headspace Vials	NKSVB-40-20	100pcs/bag
	20-60ml Vials Rack 40 Wells, Well Dia 28mm, for ND24 Screw Vials	NKSVB-40-60	100pcs/bag
	1.5ml/2ml Vials Rack 50 Wells, Well Dia 12mm	NKSVB-50-2	100pcs/bag
	4ml Vials Rack 50 Wells, Well Dia 16mm	NKSVB-50-4	100pcs/bag
	1.5ml/2ml Vials Rack 50 Wells (Thickened Type), Well Dia 12mm	NKVR-2-T	100pcs/bag

Manual / Electronic Crimpers and Decappers (11 mm / 20 mm)

Precise Fit & High Compatibility:

Fits 11mm/20mm standard vials, screw crimp vials and chromatography-specific crimp vials. Excellent sealing/de-capping adhesion to avoid vial damage and seal failure.




Two Models for Diverse Needs: Manual:
Lightweight, power-free, cost-effective for small-batch/on-site/mobile use.



Description	Spec	Cat. No.	Pack
Crimper	11mm	NKCR11	1pcs/box
	13mm	NKCR13	1pcs/box
	20mm	NKCR20	1pcs/box
Decapper	11mm	NKDC11	1pcs/box
	13mm	NKDC13	1pcs/box
	20mm	NKDC20	1pcs/box
Electronic Crimpers	11mm	NKECR11	1pcs/box
	20mm	NKECR20	1pcs/box
Electronic Decappers	11mm	NKEDC11	1pcs/box
	20mm	NKEDC20	1pcs/box


Syringe Filters and Membranes

Disposable syringe filter is a fast, convenient, and reliable filtering tool commonly used in laboratories, with a beautiful, light weight, and high cleanliness. It is mainly used for sample pre-filtration, clarification, particle removal, as well as liquid and gas sterilization and filtration. It is the preferred method for filtering small amounts of HPLC and GC samples.



Syringe filter structure form:

- Upper shell
- Membrane
- Lower shell
- Plastic particle injection molding or directwelding
- Membrane material:
 - PES, MCE, PVDF, PTFE, NY, GF, PP...
- Membrane hole size:
 - 0.1, 0.22, 0.45, 0.65, 0.8, 1.0, 1.2, 3, 5μm...
- Syringe filter diameter: 4, 13, 17, 25, 30, 33, 50mm



Features of Syringe filter:

- High flow rate
- High throughput → not easily clogged
- High retention rate (>99%)
- High voltage withstand
- Small retention volume
- Various pore sizes, diameters, and materials
- Low Dissolved Matter

Filter Items				Recommend membrane	
Air filtration				hydrophobic PTFE, PP	
Liquid filtration	Chemistry and Analysis Laboratory	Aqueous solution filtration	PH 1-14	PES	
			PH 2-13	NY, PES	
			PH 4-10	PES, MCE	
	Biology Lab	General organic solvent filtration	Strong acid, strong alkali solution, strong polarity, and corrosive solvent	PH 4-10	NY, hydrophobic PTFE; hydrophobic PVDF second choice
					hydrophobic hydrophilic PTFE
		Cell culture medium	PES, MCE, CA, hydrophilic PVDF		
		Biological solutions, proteins, peptides, etc	(Low protein adsorption) PES; second hydrophilic PVDF		
environmental laboratory	Bacterial filtration detection in solution (surface water, Bottled water, etc.)	PH 4-10	MCE, CN, CA		
General Laboratory	Various pre filtration		GF, PP		

Nylon (NY)



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,Nylon	13mm	0.22um	NKSF-N-1322	100pcs/box
		0.45um	NKSF-N-1345	100pcs/box
	25mm	1um	NKSF-N-2501	100pcs/box
		0.22um	NKSF-N-2522	100pcs/box
		0.45um	NKSF-N-2545	100pcs/box
	30mm	1um	NKSF-N-3001	100pcs/box
		0.22um	NKSF-N-3022	100pcs/box
		0.45um	NKSF-N-3045	100pcs/box

Polyethersulfone (PES)



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,PES	13mm	0.22um	NKSF-P-1322	100pcs/box
		0.45um	NKSF-P-1345	100pcs/box
	25mm	1um	NKSF-P-2501	100pcs/box
		0.22um	NKSF-P-2522	100pcs/box
		0.45um	NKSF-P-2545	100pcs/box
	30mm	1um	NKSF-P-3001	100pcs/box
		0.22um	NKSF-P-3022	100pcs/box
		0.45um	NKSF-P-3045	100pcs/box

Mixed Cellulose Esters (MCE)



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,MCE	13mm	0.22um	NKSF-M-1322	100pcs/box
		0.45um	NKSF-M-1345	100pcs/box
	25mm	1um	NKSF-M-2501	100pcs/box
		0.22um	NKSF-M-2522	100pcs/box
		0.45um	NKSF-M-2545	100pcs/box
	30mm	1um	NKSF-M-3001	100pcs/box
		0.22um	NKSF-M-3022	100pcs/box
		0.45um	NKSF-M-3045	100pcs/box

PTFE (Hydrophilic/Hydrophobic)



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia, non sterile, Hydrophobic PTFE	13mm	0.22um	NKSF-PT-1322S	100pcs/box
		0.45um	NKSF-PT-1345S	100pcs/box
	25mm	1um	NKSF-PT-2501S	100pcs/box
		0.22um	NKSF-PT-2522S	100pcs/box
		0.45um	NKSF-PT-2545S	100pcs/box
	30mm	1um	NKSF-PT-3001S	100pcs/box
		0.22um	NKSF-PT-3022S	100pcs/box
		0.45um	NKSF-PT-3045S	100pcs/box
	Syringe filter, dia, non sterile, Hydrophilic PTFE	13mm	0.22um	NKSF-PT-1322Q
0.45um			NKSF-PT-1345Q	100pcs/box
25mm		1um	NKSF-PT-2501Q	100pcs/box
		0.22um	NKSF-PT-2522Q	100pcs/box
		0.45um	NKSF-PT-2545Q	100pcs/box
30mm		1um	NKSF-PT-3001Q	100pcs/box
		0.22um	NKSF-PT-3022Q	100pcs/box
		0.45um	NKSF-PT-3045Q	100pcs/box

PVDF (Hydrophilic/Hydrophobic)



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,Hydrophobic PVDF	13mm	0.22um	NKSF-PV-1322S	100pcs/box
		0.45um	NKSF-PV-1345S	100pcs/box
	25mm	1um	NKSF-PV-2501S	100pcs/box
		0.22um	NKSF-PV-2522S	100pcs/box
		0.45um	NKSF-PV-2545S	100pcs/box
	30mm	1um	NKSF-PV-3001S	100pcs/box
		0.22um	NKSF-PV-3022S	100pcs/box
		0.45um	NKSF-PV-3045S	100pcs/box
	Syringe filter, dia , non sterile, Hydrophilic PVDF	13mm	0.22um	NKSF-PV-1322Q
0.45um			NKSF-PV-1345Q	100pcs/box
25mm		1um	NKSF-PV-2501Q	100pcs/box
		0.22um	NKSF-PV-2522Q	100pcs/box
		0.45um	NKSF-PV-2545Q	100pcs/box
30mm		1um	NKSF-PV-3001Q	100pcs/box
		0.22um	NKSF-PV-3022Q	100pcs/box
		0.45um	NKSF-PV-3045Q	100pcs/box

CA



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,CA	13mm	0.22um	NKSF-CA-1322	100pcs/box
		0.45um	NKSF-CA-1345	100pcs/box
	25mm	1um	NKSF-CA-2501	100pcs/box
		0.22um	NKSF-CA-2522	100pcs/box
		0.45um	NKSF-CA-2545	100pcs/box

CN



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,CN	13mm	0.22um	NKSF-CN-1322	100pcs/box
		0.45um	NKSF-CN-1345	100pcs/box
	25mm	1um	NKSF-CN-2501	100pcs/box
		0.22um	NKSF-CN-2522	100pcs/box
		0.45um	NKSF-CN-2545	100pcs/box

PP



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,PP	13mm	0.22um	NKSF-PP-1322	100pcs/box
		0.45um	NKSF-PP-1345	100pcs/box
	25mm	1um	NKSF-PP-2501	100pcs/box
		0.22um	NKSF-PP-2522	100pcs/box
		0.45um	NKSF-PP-2545	100pcs/box

Glass fiber



Description	diameter	Pore size	Cat. No.	Pack
Syringe filter, dia , non sterile,Glass fiber	13mm	0.22um	NKSF-GF-1322	100pcs/box
		0.45um	NKSF-GF-1345	100pcs/box
	25mm	1um	NKSF-GF-2501	100pcs/box
		0.22um	NKSF-GF-2522	100pcs/box
		0.45um	NKSF-GF-2545	100pcs/box

Sterilized Syringe Filters


To sterilize various small volume biological liquids, such as proteins, culture media, additives, buffers, reagents or drugs. 0.22um filter is used to remove all bacteria, 0.45um filter is used to remove particles from the liquid.

- Single-use, γ Radiation sterilization
- Sterilization and filtration efficiency (0.22 μ m) LRV > 7 Shortwave Monomonas,
- ATCC19146, ASTM Bacterial Challenge Experiment
- Dnase Rnase free, Non-pyrogenic, no surfactant, extremely low extract
- Extremely low protein adsorption rate < 1.5% (PES)
- Ultra high flow rate and flux, faster and more liquid filtration



Description	Specification	Part no.	packs/ctn
syringe filter, MCE, sterile	diameter 13mm, pore size 0.22um	NKSF-M-1322-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-M-1345-S-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.22um	NKSF-M-2522-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-M-2545-S	1pcs/pk, 100pk/box
syringe filter, PES, sterile	diameter 13mm, pore size 0.22um	NKSF-P-1322-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-P-1345-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.22um	NKSF-P-2522-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-P-2545-S	1pcs/pk, 100pk/box
syringe filter, Nylon, sterile	diameter 13mm, pore size 0.22um	NKSF-N-1322-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-N-1345-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.22um	NKSF-N-2522-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-N-2545-S	1pcs/pk, 100pk/box
syringe filter, hydrophobic PTFE, sterile	diameter 13mm, pore size 0.22um	NKSF-PT-1322S-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-PT-1345S-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.22um	NKSF-PT-2522S-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PT-2545S-S	1pcs/pk, 100pk/box
syringe filter, hydrophilic PTFE, sterile	diameter 13mm, pore size 0.22um	NKSF-PT-1322Q-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-PT-1345Q-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PT-2522Q-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PT-2545Q-S	1pcs/pk, 100pk/box
syringe filter, hydrophobic PVDF, sterile	diameter 13mm, pore size 0.22um	NKSF-PV-1322S-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-PV-1345S-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.22um	NKSF-PV-2522S-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PV-2545S-S	1pcs/pk, 100pk/box
syringe filter, hydrophilic PVDF, sterile	diameter 13mm, pore size 0.22um	NKSF-PV-1322Q-S	1pcs/pk, 100pk/box
	diameter 13mm, pore size 0.45um	NKSF-PV-1345Q-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PV-2522Q-S	1pcs/pk, 100pk/box
	diameter 25mm, pore size 0.45um	NKSF-PV-2545Q-S	1pcs/pk, 100pk/box

Syringless Filters Vial





FEATURES & BENEFITS

- Adapted equipment: Waters, Agilent, Shimadzu, and other mainstream chromatography instruments. Standard size suitable for auto-sampler
- Rapid sample preparation, single step process, filtering with a plunger in the vial
- Excellent cleanliness for better chromatography
- Preslit cap ensures easy and clean sample transfer
- Replace syringe, syringe filter, glass vial and cap, cost effective
- Color-coded for your convenience
- High filtration efficiency >98

SPECIFICATION

	Housing	PP
	Septa	Red PTFE/white Silicon
	Size	φ12.3*34mm
	Sterile	No
	Filling vol	480ul
	Dead vol	30ul
	Pressure	2.5 bar (36 psi)
	Filtration Area	0.28cm ²
	Max operation temp	121 °F (50 °C)
	Packaging	100/pk



Description	Specification	Part no.	packs/ctn
Filter vial, clear	Filter vial, Nylon membrane, 0.22μm, clear	NKFV-N-22-C	100pcs/pk
	Filter vial, Nylon membrane, 0.45μm, clear	NKFV-N-45-C	100pcs/pk
	Filter vial, PTFE membrane, 0.22μm, clear	NKFV-PT-22-C	100pcs/pk
	Filter vial, PTFE membrane, 0.45μm, clear	NKFV-PT-45-C	100pcs/pk
	Filter vial, PVDF Membrane, 0.22μm, clear	NKFV-PV-22-C	100pcs/pk
	Filter vial, PVDF Membrane, 0.45μm, clear	NKFV-PV-45-C	100pcs/pk
	Filter vial, PES membrane, 0.22μm, clear	NKFV-P-22-C	100pcs/pk
	Filter vial, PES membrane, 0.45μm, clear	NKFV-P-45-C	100pcs/pk
	Filter vial, RC membrane, 0.22μm, clear	NKFV-R-22-C	100pcs/pk
	Filter vial, RC membrane, 0.45μm, clear	NKFV-R-45-C	100pcs/pk

Syringless Filters Vial

Description	Specification	Part no.	packs/ctn
Filter vial, clear	Filter vial, Nylon membrane, 0.22µm, clear	NKFV-N-22-C	100pcs/pk
	Filter vial, Nylon membrane, 0.45µm, clear	NKFV-N-45-C	100pcs/pk
	Filter vial, PTFE membrane, 0.22µm, clear	NKFV-PT-22-C	100pcs/pk
	Filter vial, PTFE membrane, 0.45µm, clear	NKFV-PT-45-C	100pcs/pk
	Filter vial, PVDF Membrane, 0.22µm, clear	NKFV-PV-22-C	100pcs/pk
	Filter vial, PVDF Membrane, 0.45µm, clear	NKFV-PV-22-C	100pcs/pk
	Filter vial, PES membrane, 0.22µm, clear	NKFV-P-22-C	100pcs/pk
	Filter vial, PES membrane, 0.45µm, clear	NKFV-P-45-C	100pcs/pk
	Filter vial, RC membrane, 0.22µm, clear	NKFV-R-22-C	100pcs/pk
	Filter vial, RC membrane, 0.45µm, clear	NKFV-R-45-C	100pcs/pk
Filter vial, amber	Filter vial, Nylon membrane, 0.22µm, amber	NKFV-N-22-A	100pcs/pk
	Filter vial, Nylon membrane, 0.45µm, amber	NKFV-N-45-A	100pcs/pk
	Filter vial, PTFE membrane, 0.22µm, amber	NKFV-PT-22-A	100pcs/pk
	Filter vial, PTFE membrane, 0.45µm, amber	NKFV-PT-45-A	100pcs/pk
	Filter vial, PVDF Membrane, 0.22µm, amber	NKFV-PV-22-A	100pcs/pk
	Filter vial, PVDF Membrane, 0.45µm, amber	NKFV-PV-22-A	100pcs/pk
	Filter vial, PES membrane, 0.22µm, amber	NKFV-P-22-A	100pcs/pk
	Filter vial, PES membrane, 0.45µm, amber	NKFV-P-45-A	100pcs/pk
	Filter vial, RC membrane, 0.22µm, amber	NKFV-R-22-A	100pcs/pk
	Filter vial, RC membrane, 0.45µm, amber	NKFV-R-45-A	100pcs/pk
Filter vial, clear, preslit	Filter vial, Nylon membrane, 0.22µm, clear	NKFV-N-22-C-P	100pcs/pk
	Filter vial, Nylon membrane, 0.45µm, clear	NKFV-N-45-C-P	100pcs/pk
	Filter vial, PTFE membrane, 0.22µm, clear	NKFV-PT-22-C-P	100pcs/pk
	Filter vial, PTFE membrane, 0.45µm, clear	NKFV-PT-45-C-P	100pcs/pk
	Filter vial, PVDF Membrane, 0.22µm, clear	NKFV-PV-22-C-P	100pcs/pk
	Filter vial, PVDF Membrane, 0.45µm, clear	NKFV-PV-22-C-P	100pcs/pk
	Filter vial, PES membrane, 0.22µm, clear	NKFV-P-22-C-P	100pcs/pk
	Filter vial, PES membrane, 0.45µm, clear	NKFV-P-45-C-P	100pcs/pk
	Filter vial, RC membrane, 0.22µm, clear	NKFV-R-22-C-P	100pcs/pk
	Filter vial, RC membrane, 0.45µm, clear	NKFV-R-45-C-P	100pcs/pk
Filter vial, amber, preslit	Filter vial, Nylon membrane, 0.22µm, amber	NKFV-N-22-A-P	100pcs/pk
	Filter vial, Nylon membrane, 0.45µm, amber	NKFV-N-45-A-P	100pcs/pk
	Filter vial, PTFE membrane, 0.22µm, amber	NKFV-PT-22-A-P	100pcs/pk
	Filter vial, PTFE membrane, 0.45µm, amber	NKFV-PT-45-A-P	100pcs/pk
	Filter vial, PVDF Membrane, 0.22µm, amber	NKFV-PV-22-A-P	100pcs/pk
	Filter vial, PVDF Membrane, 0.45µm, amber	NKFV-PV-22-A-P	100pcs/pk
	Filter vial, PES membrane, 0.22µm, amber	NKFV-P-22-A-P	100pcs/pk
	Filter vial, PES membrane, 0.45µm, amber	NKFV-P-45-A-P	100pcs/pk
	Filter vial, RC membrane, 0.22µm, amber	NKFV-R-22-A-P	100pcs/pk
	Filter vial, RC membrane, 0.45µm, amber	NKFV-R-45-A-P	100pcs/pk

Membranes (25mm / 47 mm / 50 mm)

Nylon

Suitable for filtering aqueous solutions and most organic solvents. Also use with a wide range of biological preparations and can be used where other membranes are unsuitable. It is hydrophilic, eliminating the need for wetting agents that could be extracted when filtering aqueous solutions. It is flexible, durable and tear resistant.

PES(Polyethersulfone)

It is hydrophilic. Low protein binding and stable in alkaline pH. Good for aqueous applications and biological samples. It is also an alternative choice which a filter that exhibits extremely low protein binding characteristics.

PTFE

PTFE is the membrane can use with the aggressive solvents, liquids and gases that can attack other membranes. It is resistant to most acids, alkalis and solvents. The hydrophobic PTFE has applications for air and gas sterilization. Hydrophilic PTFE can be used to filter aqueous and organic mixed liquids.

MCE (Mixed Cellulose Ester)

They are modified to offer the lowest binding filters. Due to the extreme low binding characteristics, these filters provide higher throughputs than competitive offerings and reduce filter changes when filtering.

Hydrophilic PVDF

Suitable for such as multi-plexed bead-based assays, dilute protein samples, cold and hot water for injection (WFI), chemicals, aggressive solvents and sanitizing agents. Not suitable for acetone, DMSO or based >6N


	C=Compatibility	LC=Limited Compatibility	IC=Incompatibility	ND=No Data at present	
	Chemical				PVDF PES
ACIDS	acetic acid 25%		C		C C
	hydrochloric acid 25%		IC		C C
	Sulfuric acid,25%		IC		C C
	Nitric acid,25%		IC		C C
	Phosphoric acid,25%		IC		ND ND
	Formic acid,25%		IC		ND ND
	trichloroacetic acid,10%		IC		ND ND
Alcohol	Methyl alcohol,98%		C		C C
	ethyl alcohol.70%		LC		C C
	Isopropyl alcoho		C		C C
	Butanol		C		C C
	benzyl alcohol		C		C ND
	glycerinum		C		C C
Acid amides	Triethanolamine		C		ND ND
	Aniline		ND		ND ND
	Pyridine		C		IC IC
	Acetonitrile		C		C LC
Lipid	Ethyl acetate		C		C IC
	N-butyl acetate		C		IC IC
	Propyl acetate		C		IC IC
	Cellosolve acetate		ND		ND IC
	Methyl Cellosolve		ND		ND IC
	Isopropyl myristate		C		ND IC
Ketone	Acetone		C		IC IC
	Cyclohexanone		C		IC IC
	Methyl ethyl ketone		C		LC IC
	Methyl-Isobutyl Ketone		ND		LC IC
Organic Peroxide	Dioxane		C		LC IC
	Tetrahydrofuran		C		LC IC
	Dimethyl sulfoxide		C		IC IC
	sopropyl ether		ND		C C
Sundry	Phenol solution,10%		ND		LC IC
	Formaldehyde Solution,30%		C		C C
	Hydrogen Peroxide		C		ND ND
	Silicone oil		ND		C C



Description	Specification	Part no.	packs/ctn
MCE	diameter 47mm, pore size 0.22um	NKFM-M-4722	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-M-4745	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-M-5022	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-M-5045	100pcs/pk
Nylon	diameter 47mm, pore size 0.22um	NKFM-N-4722	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-N-4745	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-N-5022	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-N-5045	100pcs/pk
PES	diameter 47mm, pore size 0.22um	NKFM-P-4722	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-P-4745	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-P-5022	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-P-5045	100pcs/pk
Hydrophilic PTFE	diameter 47mm, pore size 0.22um	NKFM-PT-4722Q	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-PT-4745Q	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-PT-5022Q	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-PT-5045Q	100pcs/pk
Hydrophobic PTFE	diameter 47mm, pore size 0.22um	NKFM-PT-4722S	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-PT-4745S	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-PT-5022S	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-PT-5045S	100pcs/pk
Hydrophilic PVDF	diameter 47mm, pore size 0.22um	NKFM-PV-4722Q	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-PV-4745Q	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-PV-5022Q	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-PV-5045Q	100pcs/pk
Hydrophobic PVDF	diameter 47mm, pore size 0.22um	NKFM-PV-4722S	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-PV-4745S	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-PV-5022S	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-PV-5045S	100pcs/pk
CA	diameter 47mm, pore size 0.22um	NKFM-CA-4722	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-CA-4745	100pcs/pk
CN	diameter 47mm, pore size 5um	NKFM-CN-4705	100pcs/pk
	diameter 47mm, pore size 0.22um	NKFM-CN-4722	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-CN-4745	100pcs/pk
MCE material, sterile, individual pack	diameter 47mm, pore size 0.22um	NKFM-M-4722-S	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-M-4745-S	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-M-5022-S	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-M-5045-S	100pcs/pk
CN Grided filter membrane, sterile, individual pack	diameter 47mm, pore size 0.22um	NKFM-CN-4722G-S	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-CN-4745G-S	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-CN-5022G-S	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-CN-5045G-S	100pcs/pk
MCE Grided filter membrane, sterile, individual pack	diameter 47mm, pore size 0.22um	NKFM-M-4722G-S	100pcs/pk
	diameter 47mm, pore size 0.45um	NKFM-M-4745G-S	100pcs/pk
	diameter 50mm, pore size 0.22um	NKFM-M-5022G-S	100pcs/pk
	diameter 50mm, pore size 0.45um	NKFM-M-5045G-S	100pcs/pk

SPE column

The elution strength gradually increases during the experiments. In general, weak solvent is used during activating and loading . Weak solvent (might mixed with small amount of strong solvent) used during washing. Strong solvent (might mixed with small amount of weak solvent) during elution.



Polymeric SPE: HLB, MCX, MAX, WCX, WAX, PS/DVB

Silica based SPE: C18E, C18, C8, Phenyl, CN, NH₂, PSA, SCX, SAX, WAX, PRS, Silica, Diol

Inorganic SPE: Florisil PR, GraphiCarb, Alumina-A/N/B, Celite


Mixed mode: GraphiCarb/NH₂, C8/SCX, GraphiCarb/PSA, GraphiCar/SAX, SAX/PSA

Special SPE: Bap, Sudan Red, Plasticizer, Tea leaf, AZO Dyes, Chinese Herbal

Immunoaffinity column: aflatoxin, ochratoxin, T2, zearalenone, vomitoxin, fumonisin

IC pretreatment column: IC-H, IC-Ag, IC-Ag/H, IC-Ba, IC-Na, IC-RP

SPE Cartridge Sorbent Selection Guide




Mechanism of Action

Type	Reversed Phase	Normal Phase	Ion Exchange	Other
Mechanism	The reversed phase column sorbent is less polar substance, which extract medium or less polar substances from the extract through hydrophobic interaction.	The normal phase column sorbent is a strong polar substance, which extract strong polar substance from the extract through hydrophilic interaction.	The ion exchange column sorbent is charged, which extract corresponding charged ions by interaction force.	Solid Supported Liquid,Liquid Extraction (SLE), Graphite-Carbon Black (GCB), etc.
Typical Application	Degreasing by reversed phase column in the detection of veterinary drug residues in foods of animal origin.	Remove organic acids from plant-based foods by normal phase column in pesticide residues test.	The MCX column adsorbs the target alkaline compounds such as melamine, ractopamine, etc., then rinse to remove impurities.	Extract aromatic amines in aqueous solution by SLE in azo dyes detection.

Solvent Selection:

SPE Cartridges	Category	Low Elution Strength Solvent (Weak Solvent)	High Elution Strength Solvent (Strong Solvent)
Reversed Phase Cartridges	HLB, C18, C8	Water, Acetonitrile, Methanol	Ethyl, Acetate, Ether
Normal Phase Cartridges	Florisil, ALN, Silica	n-Hexane, Cyclohexane	Acetone, Dichloromethane
Cation Exchange Cartridges	MCX, SCX, WCX, PRS	Acidic or Neutral Solution	Alkaline Solution
Anion Exchange Cartridges	MAX, SAX, WA X	Alkaline or Neutral Solution	Acid Solution

SPE Cartridge Directions



Polymer-Based Sorbents SPE

Cat. No.	Description	Specifications
NKHLB-60-3	HLB SPE cartridge	60mg/3ml,50pk
NKHLB-150-3	HLB SPE cartridge	150mg/3ml,50pk
NKHLB-200-3	HLB SPE cartridge	200mg/3ml,50pk
NKHLB-500-3	HLB SPE cartridge	500mg/3ml,50pk
NKHLB-150-6	HLB SPE cartridge	150mg/6ml,30pk
NKHLB-200-6	HLB SPE cartridge	200mg/6ml,30pk
NKHLB-500-6	HLB SPE cartridge	500mg/6ml,30pk
NKMAX-60-3	MAX SPE cartridge	60mg/3ml,50pk
NKMAX-150-3	MAX SPE cartridge	150mg/3ml,50pk
NKMAX-150-6	MAX SPE cartridge	150mg/6ml,30pk
NKMAX-200-6	MAX SPE cartridge	200mg/6ml,30pk
NKMAX-500-6	MAX SPE cartridge	500mg/6ml,30pk

Polymer-Based Sorbents SPE

Cat. No.	Description	Specifications
NKMCX-60-3	MCX SPE cartridge	60mg/3ml,50pk
NKMCX-150-6	MCX SPE cartridge	150mg/6ml,30pk
NKMCX-200-6	MCX SPE cartridge	200mg/6ml,30pk
NKMCX-500-6	MCX SPE cartridge	500mg/6ml,30pk
NKWAX-150-3	WAX SPE cartridge	150mg/3ml,50pk
NKWAX-150-6	WAX SPE cartridge	150mg/6ml,30pk

Silica-Based Sorbents SPE

Cat. No.	Description	Specifications
NKC18-60-3	C18 SPE cartridge	60mg/3ml,50pk
NKC18-200-6	C18 SPE cartridge	200mg/6ml,30pk
NKC18-500-6	C18 SPE cartridge	500mg/6ml,30pk
NKC18-1000-6	C18 SPE cartridge	1000mg/6ml,30pk
NKPSA-500-6	PSA SPE cartridge	500mg/6ml,30pk
NKPSA-1000-6	PSA SPE cartridge	1000mg/6ml,30pk
NKSilica-500-6	Silica SPE cartridge	500mg/6ml,30pk
NKSCX-500-6	SCX SPE cartridge	500mg/6ml,30pk
NKSCX-1000-6	SCX SPE cartridge	1000mg/6ml,30pk
NKSAX-150-6	SAX SPE cartridge	150mg/6ml,30pk
NKSAX-500-6	SAX SPE cartridge	500mg/6ml,30pk
NKSAX-1000-6	SAX SPE cartridge	1000mg/6ml,30pk

Adsorption Sorbents SPE

Cat. No.	Description	Specifications
NKFL-1000-6	Florisil SPE cartridge	1000mg/6ml,30pk
NKCarb-500-6	Carb SPE cartridge	500mg/6ml,30pk
NKALN-1000-6	Alumina-N SPE cartridge	1000mg/6ml,30pk
NKALN-2000-6	Alumina-N SPE cartridge	2g/6ml,30pk
NKALN-1000-12	Alumina-N SPE cartridge	1g/12ml,20pk
NKALN-2000-12	Alumina-N SPE cartridge	2g/12ml,20pk
NKALB-1000-3	Alumina-B SPE cartridge	1000mg/3ml,50pk
NKALB-1000-6	Alumina-B SPE cartridge	1000mg/6ml,30pk
NKALB-1000-12	Alumina-B SPE cartridge	1g/12ml,20pk
NKPA-500-6	PA SPE cartridge	500mg/6ml,30pk
NKPA-1000-6	PA SPE cartridge	1000mg/6ml,30pk

Mixed-Mode Sorbents SPE

Cat. No.	Description	Specifications
NKCarb/NH2-500-500-6	Carb/NH2 SPE cartridge	500mg/500mg/6ml,30pk
NKCarb/PSA-500-500-6	Carb/PSA SPE cartridge	500mg/500mg/6ml,30pk

General Use Plastic Consumables



Universal pipette tips

Universal pipette tips are made of Polypropylene (PP) and designed to suit for single and multichannel pipettes for transferring small amounts of liquid in laboratory. And we supply DNase/RNase-free pipette tip, sterile pipette tip, non-sterile pipette tips, graduated pipette tips are suitable for most brands including Gilson, Eppendorf, Thermo, BioHit and so on.

- Highly transparent
- Smooth inner wall, reduce hanging wall residue
- Universal pipette tips and Low retention pipette tips



10ul pipette tips 32mm, 0.1-10ul, clear,

Part no.	Name	Sterile	packs/ctn
NKPT-10-B	10ul pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-10-B	10ul filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-10	10ul pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-10	10ul filter pipette tips	Yes	96pcs/rack, 50racks/ctn



10ul long pipette tips 46mm, 0.1-10ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-10-XL-B	10ul long pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-10-XL-B	10ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-10-XL	10ul long pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-10-XL	10ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



20ul filter pipette tips, 51mm, 1-20ul

Part no.	Name	Sterile	packs/ctn
NKPT-FT-20-B	20ul clear filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-20	20ul clear filter pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-20H-B	20ul yellow filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-20H	20ul yellow filter pipette tips	Yes	96pcs/rack, 50racks/ctn



20ul long filter pipette tips, 59mm, 1-20ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-FT-20-XL-B	20ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-20-XL	20ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



50ul filter pipette tips, 51mm, 1-50ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-FT-50-B	50ul clear filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-50	50ul clear filter pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-50H-B	50ul yellow filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-50H	50ul yellow filter pipette tips	Yes	96pcs/rack, 50racks/ctn



50ul long filter pipette tips, 59mm, 1-50ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-FT-50-XL-B	50ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-50-XL	50ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



100ul filter pipette tips, 51mm, 1-100ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-FT-100-B	100ul clear filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-100	100ul clear filter pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-100H-B	100ul yellow filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-100H	100ul yellow filter pipette tips	Yes	96pcs/rack, 50racks/ctn



100ul long filter pipette tips, 59mm, 1-100ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-FT-100-XL-B	100ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-100-XL	100ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



200ul pipette tips 51mm, 1-200ul, yellow

Part no.	Name	Sterile	packs/ctn
NKPT-200-B	200ul clear pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-200-B	200ul clear filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-200	200ul clear pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-200	200ul clear filter pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-200H-B	200ul yellow pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-200H-B	200ul yellow filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-200H	200ul yellow pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-200H	200ul yellow filter pipette tips	Yes	96pcs/rack, 50racks/ctn



200ul long filter pipette tips 59mm, 1-200ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-200-XL-B	200ul long pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-200-XL-B	200ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-200-XL	200ul long pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-200-XL	200ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



200ul long pipette tips 89mm, 1-200ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-200-SXL-B	200ul long pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-200-SXL-B	200ul long filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-200-SXL	200ul long pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-200-SXL	200ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn

Universal pipette tips



300ul pipette tips 52mm, 1-300ul

Part no.	Name	Sterile	packs/ctn
NKPT-300-B	300ul pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-300-B	300ul filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-300	300ul pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-300	300ul filter pipette tips	Yes	96pcs/rack, 50racks/ctn



300ul long filter pipette tips 59mm, 1-300ul

Part no.	Name	Sterile	packs/ctn
NKPT-FT-300-XL-B	300ul long filter pipette tips	Yes	1000pcs/bag, 20bags/ctn
NKPT-FT-300-XL	300ul long filter pipette tips	Yes	96pcs/rack, 50racks/ctn



1000ul pipette tips 71mm, 100-1000ul

Part no.	Name	Sterile	packs/ctn
NKPT-1000-B	1000ul pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-1000-B	1000ul filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-1000	1000ul pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-FT-1000	1000ul filter pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-1000L-B	1000ul blue pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-1000L-B	1000ul filter blue pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-1000L	1000ul blue pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-FT-1000L	1000ul filter blue pipette tips	Yes	100pcs/rack, 50racks/ctn



1000ul filter pipette tips 87mm, 100-1000ul

Part no.	Name	Sterile	packs/ctn
NKPT-1000-XL-B	1000ul pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-1000-XL-B	1000ul filter pipette tips	No	1000pcs/bag, 5bags/ctn
NKPT-1000-XL	1000ul pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-1000-XL	1000ul filter pipette tips	Yes	96pcs/rack, 50racks/ctn



1000ul pipette tips 78mm, 100-1000ul

Part no.	Name	Sterile	packs/ctn
NKPT-1000-C-B	1000ul pipette tips	Yes	1000pcs/bag, 10bags/ctn
NKPT-FT-1000-C-B	1000ul filter pipette tips	Yes	1000pcs/bag, 10bags/ctn
NKPT-1000-C	1000ul pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-FT-1000-C	1000ul filter pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-1000-C1	1000ul pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-FT-1000-C1	1000ul filter pipette tips	Yes	100pcs/rack, 50racks/ctn
NKPT-1000L-C-B	1000ul blue pipette tips	Yes	1000pcs/bag, 10bags/ctn
NKPT-FT-1000L-C-B	1000ul filter blue pipette tips	No	1000pcs/bag, 10bags/ctn
NKPT-1000L-C	1000ul blue pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-1000L-C	1000ul filter blue pipette tips	Yes	96pcs/rack, 50racks/ctn



1250ul pipette tips 102mm, 100-1250ul, clear

Part no.	Name	Sterile	packs/ctn
NKPT-1250-B	1250ul pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-FT-1250-B	1250ul filter pipette tips	No	1000pcs/bag, 20bags/ctn
NKPT-1250	1250ul pipette tips	Yes	96pcs/rack, 50racks/ctn
NKPT-FT-1250	1250ul filter pipette tips	Yes	96pcs/rack, 50racks/ctn



5ml pipette tips, clear, for Eppendorf, Nokin, Dlab, Biohit pipette

Part no.	Name	Sterile	packs/ctn
NKPT-5K-E-B	5ml pipette tips	No	300pcs/bag, 10bags/ctn
NKPT-FT-5K-E-B	5ml filter pipette tips	No	300pcs/bag, 10bags/ctn
NKPT-5K-E	5ml pipette tips	Yes	28pcs/rack, 42racks/ctn
NKPT-FT-5K-E	5ml filter pipette tips	Yes	28pcs/rack, 42racks/ctn



5ml pipette tips graduated, clear, for Eppendorf, Nokin, Dlab, Biohit pipette

Part no.	Name	Sterile	packs/ctn
NKPT-5K-G-B	5ml graduated pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-FT-5K-G-B	5ml graduated filter pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-5K-G	5ml graduated pipette tips	Yes	24pcs/rack, 30racks/ctn
NKPT-FT-5K-G	5ml filter pipette tips	Yes	24pcs/rack, 30racks/ctn



5ml pipette tips, clear, for Thermofisher, Brand, Finland pipette

Part no.	Name	Sterile	packs/ctn
NKPT-5K-T-B	5ml pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-FT-5K-T-B	5ml filter pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-5K-T	5ml pipette tips	Yes	40pcs/rack, 30racks/ctn
NKPT-FT-5K-T	5ml filter pipette tips	Yes	40pcs/rack, 30racks/ctn



10ml filter pipette tips, clear, for Eppendorf, Gilson, Dlab, Biohit, NOKE

Part no.	Name	Sterile	packs/ctn
NKPT-10K-E-B	10ml pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-FT-10K-E-B	10ml filter pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-10K-E	10ml pipette tips	Yes	24pcs/rack, 30racks/ctn
NKPT-FT-10K-E	10ml filter pipette tips	Yes	24pcs/rack, 30racks/ctn



10ml pipette tips, clear, for Thermofisher, Brand

Part no.	Name	Sterile	packs/ctn
NKPT-10K-T-B	10ml pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-FT-10K-T-B	10ml filter pipette tips	No	100pcs/bag, 10bags/ctn
NKPT-10K-T	10ml pipette tips	Yes	24pcs/rack, 30racks/ctn
NKPT-FT-10K-T	10ml filter pipette tips	Yes	24pcs/rack, 30racks/ctn

NKPT-FT-1000L-C 1000ul filter blue pipette tips Yes 96pcs/rack, 50racks/ctn

Rainin tips

Imported polypropylene (PP), USP Class VI standards

Matched with Rainin full-automatic liquid transfer workstation and LTS manual liquid transfer device, mainly used for liquid distribution and transfer, applicable to the LTSplatform.

Smooth inner wall, reduce liquid residue, ensure the accuracy of liquid suction DNA-free, DNase-free, RNase-free, Pyrogen/Endotoxin free
Ultra-grinding & molding technology, release agent free

Withstand temperature ranging from -80°C to 121 °C, stable performance under high temperature and high pressure.

Extra thin-wall design increase product flexibility and improve heat transfer efficiency



Part no.	Description	packs/ctn
NKRN-20-B	RN 20ul Clear, Bulk, bag pack	1000pcs/bag, 10bags/ctn
NKRN-20	RN 20ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-20-L-B	RN 20ul Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 10bags/ctn
NKRN-20-L	RN 20ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-FT-20-B	RN 20ul Filter, Clear, Bulk, bag pack	1000pcs/bag, 10bags/ctn
NKRN-FT-20	RN 20ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-FT-20-L	RN 20ul Filter, Low Retention, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-FT-20-L-B	RN 20ul Filter, Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 10bags/ctn
NKRN-20-RE	RN 20ul, Reload, Clear	960pcs/pack, 10packs/ctn
NKRN-20-L-RE	RN 20ul, Reload, Low Retention, Clear	960pcs/pack, 10packs/ctn

Part no.	Description	packs/ctn
NKRN-200-B	RN 200ul Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-200	RN 200ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-200-L-B	RN 200ul Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-200-L	RN 200ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-FT-200-B	RN 200ul Filter, Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-FT-200	RN 200ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKRN-FT-200-L	RN 200ul Filter, Low Retention, Clear, Sterile, Rack	96pcs/rack, 40racks/ctn
NKRN-FT-200-L-B	RN 200ul Filter, Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-200-RE	RN 200ul, Reload, Clear	960pcs/pack, 10packs/ctn
NKRN-200-L-RE	RN 200ul, Reload, Low Retention, Clear	960pcs/pack, 10packs/ctn

Part no.	Description	packs/ctn
NKRN-300-B	RN 300ul Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-300	RN 300ul Clear, Sterile, Rack	1000pcs/bag, 20bags/ctn
NKRN-300-L-B	RN 300ul Low Retention, Clear, Bulk, bag pack	96pcs/rack, 40racks/ctn
NKRN-300-L	RN 300ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 40racks/ctn
NKRN-FT-300-B	RN 300ul Filter, Clear, Bulk, bag pack	1000pcs/bag, 20bags/ctn
NKRN-FT-300	RN 300ul Filter, Clear, Sterile, Rack	1000pcs/bag, 5bags/ctn
NKRN-FT-300-L	RN 300ul Filter, Low Retention, Clear, Sterile, Rack	96pcs/rack, 40racks/ctn
NKRN-FT-300-L-B	RN 300ul Filter, Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 5bags/ctn
NKRN-300-RE	RN 300ul, Reload, Clear	768pcs/pack, 10packs/ctn
NKRN-300-L-RE	RN 300ul, Reload, Low Retention, Clear	768pcs/pack, 10packs/ctn

Part no.	Description	packs/ctn
NKRN-1000-B	RN 1000ul Tips, Clear, Bulk, bag pack	96pcs/rack, 40racks/ctn
NKRN-1000	RN 1000ul Tips, Clear, Sterile, Rack	1000pcs/bag, 5bags/ctn
NKRN-1000-L-B	RN 1000ul Low Retention, Clear, Bulk, bag pack	96pcs/rack, 40racks/ctn
NKRN-1000-L	RN 1000ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 40racks/ctn
NKRN-FT-1000-B	RN 1000ul Filter, Clear, Bulk, bag pack	1000pcs/bag, 5bags/ctn
NKRN-FT-1000	RN 1000ul Tips, Filter, Clear, Sterile, Rack	1000pcs/bag, 5bags/ctn
NKRN-FT-1000-L	RN 1000ul Filter, Low Retention, Clear, Sterile, Rack	96pcs/rack, 40racks/ctn
NKRN-FT-1000-L-B	RN 1000ul Filter, Low Retention, Clear, Bulk, bag pack	1000pcs/bag, 5bags/ctn
NKRN-1000-RE	RN 1000ul, Reload, Clear	768pcs/pack, 10packs/ctn
NKRN-1000-L-RE	RN 1000ul, Reload, Low Retention, Clear	768pcs/pack, 10packs/ctn

Beckman tips



Part no.	Description	packs/ctn
NKBK-20	BKM 20ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-FT-20	BKM 20ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-50	BKM 50ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-50-L	BKM 50ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-FT-50	BKM 50ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-FT-50-L	BKM 50ul Low Retention, Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-250	BKM 250ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-250-L	BKM 250ul Low Retention, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-FT-50	BKM 250ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKBK-FT-250-L	BKM 250ul Low Retention, Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn

Hamilton tips



Part no.	Description	packs/ctn
NKHM-50B	HMT 50ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-FT-50B	HMT 50ul Filter, Conductive, Black, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-50	HMT 50ul Clear, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-FT-50	HMT 50ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-300B	HMT 300ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-FT-300B	HMT 300ul Filter, Conductive, Black, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-300	HMT 300ul Clear, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-FT-300	HMT 300ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 60racks/ctn
NKHM-1000B	HMT 1000ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 40racks/ctn
NKHM-FT-1000B	HMT 1000ul Filter, Conductive, Black, Sterile, Blister Rack	96pcs/rack, 40racks/ctn
NKHM-1000	HMT 1000ul Clear, Sterile, Blister Rack	96pcs/rack, 40racks/ctn
NKHM-FT-1000	HMT 1000ul Tips, Filter, Clear, Sterile, Blister Rack	96pcs/rack, 40racks/ctn

Tecan tips



Part no.	Description	packs/ctn
NKTC-10B	TCN 10ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-10B	TCN 10ul Conductive, Filter, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-10B-1	TCN 10ul Conductive, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-10B-1	TCN 10ul Conductive, Filter, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-50B	TCN 50ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-50B	TCN 50ul Conductive, Filter, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-50B-1	TCN 50ul Conductive, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-50B-1	TCN 50ul Conductive, Filter, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-200B	TCN 200ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-200B	TCN 200ul Conductive, Filter, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-200B-1	TCN 200ul Conductive, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-200B-1	TCN 200ul Conductive, Filter, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-1000B	TCN 1000ul Conductive, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-1000B	TCN 1000ul Conductive, Filter, Black, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-1000B-1	TCN 1000ul Conductive, Black, Sterile, Rack	96pcs/rack, 50racks/ctn































Part no.	Description	packs/ctn
NKTC-FT-1000B-1	TCN 1000ul Conductive, Filter, Black, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-50-1	TCN 50ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-50-1	TCN 50ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-50	TCN 50ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-50	TCN 50ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-200-1	TCN 200ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-200-1	TCN 200ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-200	TCN 200ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-200	TCN 200ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-1000-1	TCN 1000ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-1000-1	TCN 1000ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-1000	TCN 1000ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-1000	TCN 1000ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn

Part no.	Description	packs/ctn
NKTC-50-1	TCN 50ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-50-1	TCN 50ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-50	TCN 50ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-50	TCN 50ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-200-1	TCN 200ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-200-1	TCN 200ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-200	TCN 200ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-200	TCN 200ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-1000-1	TCN 1000ul Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-FT-1000-1	TCN 1000ul Filter, Clear, Sterile, Rack	96pcs/rack, 50racks/ctn
NKTC-1000	TCN 1000ul Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn
NKTC-FT-1000	TCN 1000ul Filter, Clear, Sterile, Blister Rack	96pcs/rack, 24racks/ctn

Centrifuge tubes

- Open and close conveniently with one hand
- Special safe-lock cap, tight-fitting lid provides a precise and safety sealing
- Thin membrane lid can be pierced by a syringe
- Accurate molded graduation for volume estimation
- Flat lid and large frosted area on tube for easy sample identification
- Guaranteed functionality between -86°C to 121°C




 <p>NKMCT-005-S centrifuge tube, 0.5ml, clear 1000pcs/bag, 25bag/ctn</p>	 <p>NKMCT-015-S centrifuge tube, 1.5ml, clear, snap cap 500pcs/bag, 20bag/ctn</p>	 <p>NKMCT-015-F centrifuge tube, 1.5ml, clear, flat cap 500pcs/bag, 20bag/ctn</p>	 <p>NKMCT-015-L centrifuge tube, 1.5ml, clear, safe lock 500pcs/bag, 10bag/ctn</p>
 <p>NKMCT-020-S centrifuge tube, 2.0ml, clear, snap cap 500pcs/bag, 20bag/ctn</p>	 <p>NKMCT-020-F centrifuge tube, 2.0ml, clear, flat cap 500pcs/bag, 20bag/ctn</p>	 <p>NKMCT-020-F-C centrifuge tube, 2.0ml, clear, flat cap, (C type) 500pcs/bag, 10bag/ctn</p>	 <p>NKMCT-020-L centrifuge tube, 2.0ml, clear, safe lock 500pcs/bag, 10bag/ctn</p>
 <p>NKCT-5-R-F centrifuge tube, 5.0ml, clear, round bottom, flat cap 300pcs/bag, 20bag/ctn</p>	 <p>NKCT-5-C-F centrifuge tube, 5.0ml, clear, conical bottom, flat cap 100pcs/bag, 20bag/ctn</p>	 <p>NKCT-5-C-S centrifuge tube, 5.0ml, clear, conical bottom, screw cap 100pcs/bag, 20bag/ctn</p>	 <p>NKCT-5-C-L centrifuge tube, 5.0ml, clear, conical bottom, safe lock cap 100pcs/bag, 20bag/ctn</p>
 <p>NKCT-10-C centrifuge tube, 10ml, conical bottom 100pcs/bag, 20bag/ctn</p>	 <p>NKCT-10-R centrifuge tube, 10ml, round bottom 100pcs/bag, 20bag/ctn</p>	 <p>NKCT-15-R centrifuge tube, 15ml, round bottom 100pcs/bag, 20bag/ctn</p>	 <p>NKCT-15-C centrifuge tube, 15ml, conical bottom 100pcs/bag, 20bag/ctn</p>
 <p>NKCT-50-C centrifuge tube, 50ml, conical bottom 50pcs/bag, 20bag/ctn</p>	 <p>NKCT-50-S centrifuge tube, 50ml, self stand bottom 50pcs/bag, 10bag/ctn</p>	 <p>NKCT-50-R centrifuge tube, 50ml, round bottom 50pcs/bag, 20bag/ctn</p>	 <p>NKMCT-015-B centrifuge tube, 1.5ml, flat cap(B type 25000xg), Dnase&Rnase free, sterile 50pcs/bag, 20bag/ctn</p>
 <p>NKMCT-020-B centrifuge tube, 2.0ml, flat cap(B type 25000xg), Dnase&Rnase free, sterile 50pcs/bag, 10bag/ctn</p>	 <p>NKCT-10-B centrifuge tube, 10ml, conical bottom(B type 18000xg), sterile 25pcs/paper box, 20rack/ctn</p>	 <p>NKCT-15-B centrifuge tube, 15ml, conical bottom(B type 18000xg), sterile 25pcs/bag, 20bag/ctn</p>	 <p>NKCT-50-B centrifuge tube, 50ml, conical bottom(B type 18000xg), sterile 25pcs/paper box, 20rack/ctn</p>
 <p>NKCT-15-C2 centrifuge tube, 15ml, conical bottom, individual pack, sterile 50pcs/bag, 20bag/ctn</p>	 <p>NKCT-15-C1 centrifuge tube, 15ml, conical bottom, paper box 50pcs/paper box, 20rack/ctn</p>	 <p>NKCT-15-C3 centrifuge tube, 15ml, conical bottom, foam rack 50pcs/paper box, 20rack/ctn</p>	
 <p>NKCT-50-C2 centrifuge tube, 50ml, conical bottom, individual pack, sterile 25pcs/bag, 20bag/ctn</p>	 <p>NKCT-50-C1 centrifuge tube, 50ml, conical bottom, paper box 25pcs/paper box, 20rack/ctn</p>	 <p>NKCT-50-C3 centrifuge tube, 50ml, conical bottom, foam rack 25pcs/paper box, 20rack/ctn</p>	

Pasteur pipette

- Made of LDPE material, Non-toxic, adapted to draw, transfer or carry for small volumeliquid.
- Optimizing process on surface tension, easy for the liquid flowing.Hightransparency, easy to observe.
- Can be bend with a certain angle, which is convenient for drawing or adding liquid in irregular or micro container.
- High elasticity, adapted to fast liquid transfer without leakage.
- Convenient and accurate for use with good repeatability.
- Heat-seal on the tip of the pipette can achieve liquid carry.
- Available in bulk or individual pack.
- Available in E.O. or Gamma radiation sterile.

Pasteur pipette



Part no.	Description	packs/ctn
NKPP-02-B	0.2ml, 65mm, bulk	1000pcs/bag, 50000pcs/ctn
NKPP-05-B	0.5ml, 85mm, bulk	500pcs/bag, 40000pcs/ctn
NKPP-1-B-1	1ml, 145mm, bulk	100pcs/bag, 5000pcs/ctn
NKPP-1-D-S-1	1ml, 145mm, Ind. peel wrap, sterile	500pcs/box, 3000pcs/ctn
NKPP-1-B-2	1ml, 160mm, bulk	100pcs/bag, 5000pcs/ctn
NKPP-1-D-S-2	1ml, 160mm, Ind. peel wrap, sterile	500pcs/box, 3000pcs/ctn
NKPP-2-B	2ml, 150mm, bulk	100pcs/bag, 5000pcs/ctn
NKPP-2-D-S	2ml, 150mm, Ind. peel wrap, sterile	500pcs/box, 3000pcs/ctn
NKPP-3-B-1	3ml, 160mm, bulk, A type	100pcs/bag, 5000pcs/ctn
NKPP-3-D-S-1	3ml, 160mm, Ind. peel wrap, sterile, A type	500pcs/box, 3000pcs/ctn
NKPP-3-B-2	3ml, 160mm, bulk, B type	100pcs/bag, 5000pcs/ctn
NKPP-3-D-S-2	3ml, 160mm, Ind. peel wrap, sterile, B type	500pcs/box, 3000pcs/ctn
NKPP-3-B-3	3ml, 178mm, bulk	100pcs/bag, 5000pcs/ctn
NKPP-3-D-S-3	3ml, 178mm, Ind. peel wrap, sterile	200pcs/bag, 2000pcs/ctn
NKPP-5-B	5ml, 210mm, bulk	250pcs/box, 2500pcs/ctn
NKPP-5-D-S	5ml, 210mm, Ind. peel wrap, sterile	500pcs/box, 2000pcs/ctn
NKPP-10-B	10ml, 296mm, bulk	200pcs/box, 1000pcs/ctn

Serological pipette

- 1.0ml, 2.0ml, 5.0ml, 10ml pipettes are shaped by stretching method.
- 25.0ml, 50.0ml, 100ml pipettes are welded to the tube by supersonic wave.
- Different colors correspond to different specification, easy to distinguish.
- Maximally reduce the liquid adhesion to the tube inner wall and improve the sampling accuracy.
- Excellent double-side graduation design ideally meet with the demand for liquid adding or reducing.
- All serological pipettes are assembled with high quality filters.
- Strict integrity test on every pipette.
- Sterilized by radiation/E.O., non-pyrogenic and DNase/RNase-free.

Part no.	Description	packs/ctn
NKSP-1	serological pipette, 1ml, 1/100, 270mm, yellow mark, sterile, individual peel pack	50pcs/bag, 20bags/ctn
NKSP-2	serological pipette, 2ml, 1/100, 277mm, green mark, sterile, individual peel pack	50pcs/bag, 20bags/ctn
NKSP-5	serological pipette, 5ml, 1/10, 345mm, blue mark, sterile, individual peel pack	50pcs/bag, 4bags/ctn
NKSP-10	serological pipette, 10ml, 1/10, 341mm, orange mark, sterile, individual peel pack	50pcs/bag, 4bags/ctn
NKSP-25	serological pipette, 25ml, 2/10, 345mm, red mark, sterile, individual peel pack	50pcs/bag, 4bags/ctn
NKSP-50	serological pipette, 50ml, 5/10, 345mm, purple mark, sterile, individual peel pack	25pcs/bag, 4bags/ctn
NKSP-100	serological pipette, 100ml, 1, 345mm, pink mark, sterile, individual peel pack	50pcs/bag, 6bags/ctn



Test tube

The test tube with screw cap is made of PS material with smooth and smooth inner wall, excellent transparency and strong chemical resistance. PP material test tube with screw cap can withstand the centrifugal rate of 5000 revolutions per minute, no pipe wall infiltration or pipe body damage.

There are various specifications of the products to meet the requirements of different experiments.



Part no.	Description	packs/ctn
NKTT-ST-12	plastic test tube, Urinary sediment tube, 12ml	2000pcs/ctn
NKTT-1	plastic test tube, PS, 5ml, 12*75mm	5000pcs/ctn
NKTT-2	plastic test tube, PS, 10ml, 16*100mm	2000pcs/ctn
NKTT-1-CA	plastic test tube with cap, PS, 12*75mm(assembled)	500pcs/bag
NKTT-2-CA	plastic test tube with cap, PS, 16*100mm(assembled)	250pcs/bag
NKTT-3-CA	plastic test tube with cap, PS, 12*100mm(assembled)	400pcs/bag
NKTT-4-CA	plastic test tube with cap, PS, 16*150mm(assembled)	100pcs/bag
NKTT-4-CA-D-S	plastic test tube with cap, PS, 16*150mm(assembled), Ind wrap, sterile	2000pcs/ctn

Reagent bottle

- Professional leak-proof bottle mouth design, no need for inner cap or inner gasket protection, excellent sealing performance, 100% guaranteed no leakage, wide mouth design, easy access to liquids.
- Imported high-quality raw materials, excellent physical and chemical indicators, with strong resistance to pressure, impact and acid resistance, polypropylene high temperature resistance of 121 °C, high density polyethylene HDPE low temperature resistance -40 °C, and Meeting the requirements of autoclave processing.



Reagent bottle, PP, wide mouth, clear color

Part no.	Description	packs/ctn
NKRB-0008-P-W-C	8ml	150pcs/bag, 10bags/ctn
NKRB-0015-P-W-C	15ml	120pcs/bag, 10bags/ctn
NKRB-0030-P-W-C	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-P-W-C	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-P-W-C	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-P-W-C	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-P-W-C	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-P-W-C	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, PP, wide mouth, brown color

Part no.	Description	packs/ctn
NKRB-0008-P-W-B	8ml	150pcs/bag, 10bags/ctn
NKRB-0015-P-W-B	15ml	120pcs/bag, 10bags/ctn
NKRB-0030-P-W-B	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-P-W-B	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-P-W-B	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-P-W-B	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-P-W-B	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-P-W-B	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, HDPE, wide mouth, clear color

Part no.	Description	packs/ctn
NKRB-0008-H-W-C	8ml	150pcs/bag, 10bags/ctn
NKRB-0015-H-W-C	15ml	120pcs/bag, 10bags/ctn
NKRB-0030-H-W-C	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-H-W-C	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-H-W-C	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-H-W-C	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-H-W-C	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-H-W-C	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, HDPE, wide mouth, brown color

Part no.	Description	packs/ctn
NKRB-0008-H-W-B	8ml	150pcs/bag, 10bags/ctn
NKRB-0015-H-W-B	15ml	120pcs/bag, 10bags/ctn
NKRB-0030-H-W-B	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-H-W-B	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-H-W-B	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-H-W-B	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-H-W-B	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-H-W-B	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, PP, narrow mouth, clear color

Part no.	Description	packs/ctn
NKRB-0030-P-N-C	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-P-N-C	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-P-N-C	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-P-N-C	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-P-N-C	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-P-N-C	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, PP, narrow mouth, brown color

Part no.	Description	packs/ctn
NKRB-0030-P-N-B	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-P-N-B	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-P-N-B	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-P-N-B	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-P-N-B	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-P-N-B	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, HDPE, narrow mouth, clear color

Part no.	Description	packs/ctn
NKRB-0030-H-N-C	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-H-N-C	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-H-N-C	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-H-N-C	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-H-N-C	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-H-N-C	1000ml	6pcs/bag, 14bags/ctn

Reagent bottle, HDPE, narrow mouth, brown color

Part no.	Description	packs/ctn
NKRB-0030-H-N-B	30ml	100pcs/bag, 10bags/ctn
NKRB-0060-H-N-B	60ml	100pcs/bag, 10bags/ctn
NKRB-0125-H-N-B	125ml	50pcs/bag, 5bags/ctn
NKRB-0250-H-N-B	250ml	25pcs/bag, 8bags/ctn
NKRB-0500-H-N-B	500ml	12pcs/bag, 9bags/ctn
NKRB-1000-H-N-B	1000ml	6pcs/bag, 14bags/ctn

Deep well plate


English Translation (Professional & Concise for Product Specs)

Durability: Supports storage at $-40\text{ }^{\circ}\text{C} \sim -80\text{ }^{\circ}\text{C}$, withstands a maximum centrifugal force of 4000g, compatible with autoclaving at 20Psi/121 $^{\circ}\text{C}$ for 20 minutes, with uniform heating.

Material & Compatibility: Imported PP material, high stability, non-reactive with test reagents/DMSO/water, low residual liquid and heavy metal content.

Cleanliness & Safety: DNase/RNase-free and pyrogen-free certified, manufactured and packaged in a Class 100,000 cleanroom.


Usability: Three leak-proof sealing options (adhesive seals, mat caps, heat seals), excellent flatness enabling long-distance transportation, compatible with automated instruments.



Part no.	Description	packs/ctn
NKDP-96220-R-F	96 well Deep well plate, 2.2ml, round hole, F bottom	5pcs/pack, 50pcs/ctn
NKDP-96220-R-U	96 well Deep well plate, 2.2ml, round hole, U bottom	5pcs/pack, 50pcs/ctn
NKDP-96220-R-V	96 well Deep well plate, 2.2ml, round hole, V bottom	5pcs/pack, 50pcs/ctn
NKDP-96220-S-F	96 well Deep well plate, 2.2ml, square hole, F bottom	5pcs/pack, 50pcs/ctn
NKDP-96220-S-U	96 well Deep well plate, 2.2ml, square hole, U bottom	5pcs/pack, 50pcs/ctn
NKDP-96220-S-V	96 well Deep well plate, 2.2ml, square hole, V bottom	5pcs/pack, 50pcs/ctn
























PCR tube/PCR plate

Made of clean PP material, thin wall and uniform thickness can keep good heat transfer. 0.2ml strips is available in two types: Transparent and white, respectively applicable to normal PCR and real-time PCR (qPCR) reaction. Good sealing performance of the 8-strip while lid closed. Easy to open the cap without contamination. Sterilized by Gamma radiation.

 <p>NKPCR-2 PCR single tube, 0.2ml, flat cap, clear, thin wall 1000pcs/bag, 10bags/ctn</p>	 <p>NKPCR-2G PCR single tube, 0.2ml, convex cap, clear, thin wall (economical type) 1000pcs/bag, 20bags/ctn</p>	 <p>NKPCR-5 PCR single tube, 0.5ml, flat cap, clear 500pcs/bag, 10box/ctn</p>	 <p>NKPCR-1-4-C PCR 4 strips tube, 0.1ml, separate cap 250strips/bag, 10bags/ctn</p>
 <p>NKPCR-1-8-C PCR 8 strips tube, 0.1ml, clear, separate flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-1-8-W PCR 8 strips tube, 0.1ml, white, separate flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-1-8-C1 PCR 8 strips tube, 0.1ml, clear, hinged flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-1-8-C2 PCR 8 strips tube, 0.1ml, clear, separate convex cap 125strips/box, 10box/ctn</p>
 <p>NKPCR-2-8-C PCR 8 strips tube, 0.2ml, clear, separate flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-2-8-W PCR 8 strips tube, 0.2ml, white, separate flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-2-8-C1 PCR 8 strips tube, 0.2ml, clear, hinged flat cap 125strips/box, 10box/ctn</p>	 <p>NKPCR-2-8-C2 PCR 8 strips tube, 0.2ml, clear, separate convex cap 125strips/box, 10box/ctn</p>
 <p>NKPCR-P-1-N-C-P PCR plate 96 well, 0.1ml, non-skirted, clear, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-1-N-W-P PCR plate 96 well, 0.1ml, non-skirted, white, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-1-F-C-P PCR plate 96 well, 0.1ml, full-skirted, clear, printed 10pcs/box, 50pcs/ctn</p>	
 <p>NKPCR-P-1-F-W-P PCR plate 96 well, 0.1ml, full-skirted, white, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-N-C-P PCR plate 96 well, 0.2ml, non-skirted, clear, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-N-W-P PCR plate 96 well, 0.2ml, non-skirted, white, printed 10pcs/box, 50pcs/ctn</p>	
 <p>NKPCR-P-2-S-C-P PCR plate 96 well, 0.2ml, semi-skirted, clear, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-S-W-P PCR plate 96 well, 0.2ml, semi-skirted, white, printed 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-N-C PCR plate 96 well, 0.2ml, non-skirted, clear 10pcs/box, 50pcs/ctn</p>	
 <p>NKPCR-P-2-S-C PCR plate 96 well, 0.2ml, semi-skirted, clear 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-1-S-C-A PCR plate 96 well, 0.1ml, semi-skirted, clear, for ABI 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-S-C-A PCR plate 96 well, 0.2ml, semi-skirted, clear, for ABI 10pcs/box, 50pcs/ctn</p>	
 <p>NKPCR-P-1-S-C-R PCR plate 96 well, 0.1ml, semi-skirted, clear, for Roche 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-P-2-S-W-R PCR plate 96 well, 0.1ml, semi-skirted, white, for Roche 10pcs/box, 50pcs/ctn</p>	 <p>NKPCR-384 PCR plate 384 well, clear 150pcs/ctn</p>	
 <p>NKPCR-384 PCR plate 384 well, white 150pcs/ctn</p>	 <p>NKPCR-SF PCR plate sealing film, clear 100pcs/box, 40box/ctn</p>	 <p>NKPCR-R PCR tube rack, 96 well 270pcs/ctn</p>	

Cryogenic tube

- Cryogenic tube is made of medical grade polypropylene (PP) and is endotoxin-free. They are disposable laboratory consumables dedicated to storing biological samples.
- The skirt lines on the cap and the tube body are easy for screw. The cap has a built-in silicone O- ring to prevent liquid leakage.
- When operating according to the standard, there is no bursting, good sealing and no leakage.
- The tightness can be guaranteed even at the lowest storage temperature, ensuring the safety of samples and laboratories.
- Can withstand low temperature. (be stored in the gas phase liquid nitrogen and can withstand low temperature to -196°C .)
- Multiple color options for the cap.

 <p>NKC-2 cryovial, 2ml economical type, with O-ring, Temp. stand-86~121 C 500pcs/bag, 20bag/ctn</p>	 <p>NKC-5 cryovial, 5ml economical type, with O-ring, Temp. stand-86~121 C 200pcs/bag, 20bag/ctn</p>	 <p>NKC-10 cryovial, 10ml economical type, with O-ring, Temp. stand-86~121 C, cap&tube unassembled 100pcs/bag, 20bag/ctn</p>	 <p>NKC-05-B cryovial, 0.5ml, black O-ring, Temp. stand-86~121 C, cap&tube unassembled 500pcs/bag, 20bag/ctn</p>
 <p>NKC-15-B cryovial, 1.5ml, black O-ring, Temp. stand-86~121 C, cap&tube unassembled 500pcs/bag, 20bag/ctn</p>	 <p>NKC-2-B cryovial, 2.0ml, black O-ring, Temp. stand-86~121 C, cap&tube unassembled 500pcs/bag, 20bag/ctn</p>	 <p>NKC-2-I-1 cryovial, 2ml, *bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-2-E-1 cryovial, 2ml, *bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>
 <p>NKC-2-I-2 cryovial, 2ml, round bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-2-E-2 cryovial, 2ml, round bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-2-I-3 cryovial, 2ml, () bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-2-E-3 cryovial, 2ml, () bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>
 <p>NKC-3-I-1 cryovial, 3ml, *bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-3-E-1 cryovial, 3ml, *bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-4-I-1 cryovial, 4ml, *bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-4-E-1 cryovial, 4ml, *bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>
 <p>NKC-4-I-2 cryovial, 4ml, round bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-5-I-1 cryovial, 5ml, *bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-5-E-1 cryovial, 5ml, *bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-5-I-2 cryovial, 5ml, round bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>
 <p>NKC-5-E-2 cryovial, 5ml, round bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-5-I-3 cryovial, 5ml, () bottom, internal screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	 <p>NKC-5-E-3 cryovial, 5ml, () bottom, external screw cap, Temp. stand-196~121 C 100pcs/bag, 20bag/ctn</p>	

Cell culture dish

- Using pure polystyrene raw materials in line with USP-Class VI standard, GMP 100,000 clean workshop production.
- There are bumps on the inside of the upper cover to facilitate gas exchange during culture.
- The upper cover and lower bottom have corresponding stacking ring design to ensure stability when multiple products are stacked.
- The lower bottom has a serrated ring design, which makes grasping more convenient and safe, and is not easy to slip off.
- There are coordinate numbers at the bottom for easy marking and positioning
- No pyrogen, no endotoxin, no cytotoxicity.
- Irradiation Sterilization, SAL 10.⁻⁶



TC treated / non TC



Glass bottom,TC

Part no.	Description	packs/ctn
NKCD-35	cell culture dish, 35mm, non TC treated, sterile	10pcs/bag, 50bags/ctn
NKCD-60	cell culture dish, 60mm, non TC treated, sterile	10pcs/bag, 50bags/ctn
NKCD-90	cell culture dish, 90mm, non TC treated, sterile (100mm)	10pcs/bag, 30bags/ctn
NKCD-150	cell culture dish, 150mm, non TC treated, sterile	5pcs/bag, 20bags/ctn
NKCD-60-TC	cell culture dish, 60mm, TC treated, sterile	10pcs/bag, 50bags/ctn
NKCD-90-TC	cell culture dish, 90mm, TC treated, sterile (100mm)	10pcs/bag, 30bags/ctn
NKCD-150-TC	cell culture dish, 150mm, TC treated, sterile	5pcs/bag, 20bags/ctn
NKCD-35-G15-TC	cell culture dish, 35mm, glass bottom, 15mm pore size, TC treated, sterile	10pcs/bag, 50bags/ctn
NKCD-35-G20-TC	cell culture dish, 35mm, glass bottom, 20mm pore size, TC treated, sterile	10pcs/bag, 50bags/ctn

Cell culture plate

Product Feature:

- USP Class-6 Polystyrene
- Anti-slip area on base sides is easy to grip.
- Numerical codes are for easy identification.
- DNase-free, RNase-free, Non-pyrogenic, non-endotoxin, non-cytotoxicity
- Irradiation sterilized, SAL-10⁻⁶




6well

Part no.	Description	packs/ctn
NKCP-6	non TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-6-TC	TC treated, sterile	1pcs/bag,50pcs/ctn



48well

Part no.	Description	packs/ctn
NKCP-48	non TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-48-TC	TC treated, sterile	1pcs/bag,50pcs/ctn



12well

Part no.	Description	packs/ctn
NKCP-12	non TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-12-TC	TC treated, sterile	1pcs/bag,50pcs/ctn



24well

Part no.	Description	packs/ctn
NKCP-24	non TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-24-TC	TC treated, sterile	1pcs/bag,50pcs/ctn




96well


Part no.	Description	packs/ctn
NKCP-96	non TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-96-TC	TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-96-U	U bottom, nonTC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-96-U-TC	U bottom, TC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-96-V	V bottom, nonTC treated, sterile	1pcs/bag,50pcs/ctn
NKCP-96-V-TC	V bottom, TC treated, sterile	1pcs/bag,50pcs/ctn

Cell culture flask

- By pure polystyrene raw material in line with USP-Class VI standard, GMP 100,000 clean workshop production.
- 2 types of lids, sealed lids and vented lids There is a frosted writing area on the side of the culture bottle. neck and a volume scale on the side for easy identification.
- Stackable design, easy to stack, not easy to slip.
- No pyrogen, no endotoxin, no cytotoxicity.
- Irradiation Sterilization, SAL 10⁶




- ① vented cap
- ② sealed cap



Part no.	Description	packs/ctn
NKCF-25-S	cell culture flask, 25 cm ² culture area, sealed cap, non TC, sterile	10pcs/bag, 20bags/ctn
NKCF-25-V	cell culture flask, 25 cm ² culture area, vented cap, non TC, sterile	10pcs/bag, 20bags/ctn
NKCF-25-S-TC	cell culture flask, 25 cm ² culture area, sealed cap, TC treated, sterile	10pcs/bag, 20bags/ctn
NKCF-25-V-TC	cell culture flask, 25 cm ² culture area, vented cap, TC treated, sterile	10pcs/bag, 20bags/ctn
NKCF-75-S	cell culture flask, 75 cm ² culture area, sealed cap, non TC, sterile	5pcs/bag, 20bags/ctn
NKCF-75-V	cell culture flask, 75 cm ² culture area, vented cap, non TC, sterile	5pcs/bag, 20bags/ctn
NKCF-75-S-TC	cell culture flask, 75 cm ² culture area, sealed cap, TC treated, sterile	5pcs/bag, 20bags/ctn
NKCF-75-V-TC	cell culture flask, 75 cm ² culture area, vented cap, TC treated, sterile	5pcs/bag, 20bags/ctn
NKCF-175-S	cell culture flask, 175 cm ² culture area, sealed cap, non TC, sterile	5pcs/bag, 8bags/ctn
NKCF-175-V	cell culture flask, 175 cm ² culture area, vented cap, non TC, sterile	5pcs/bag, 8bags/ctn
NKCF-175-S-TC	cell culture flask, 175 cm ² culture area, sealed cap, TC treated, sterile	5pcs/bag, 8bags/ctn
NKCF-175-V-TC	cell culture flask, 175 cm ² culture area, vented cap, TC treated, sterile	5pcs/bag, 8bags/ctn
NKCF-225-S	cell culture flask, 225 cm ² culture area, sealed cap, non TC, sterile	5pcs/bag, 5bags/ctn
NKCF-225-V	cell culture flask, 225 cm ² culture area, vented cap, non TC, sterile	5pcs/bag, 5bags/ctn
NKCF-225-S-TC	cell culture flask, 225 cm ² culture area, sealed cap, TC treated, sterile	5pcs/bag, 5bags/ctn
NKCF-225-V-TC	cell culture flask, 225 cm ² culture area, vented cap, TC treated, sterile	5pcs/bag, 5bags/ctn

Cell culture flask

- The bottle is seamless and blow-moulded in one go to avoid surface knots or rupture of cells against the wall.
- The surface finish meets the requirement of mirror surface and the permeability is close to glass, reducing the risk of cell wall adhesion.
- Bottle mouth size meets international SBS standards, suitable for quick connection to bioreactors.
- Sterilised, no DNase/RNase, no pyrogens and no endotoxins. Individually vacuum packed.
- Equipped with a 0.2um gas permeable cap for better mixing of cells and culture medium, increasing air permeability and allowing cells to grow to higher densities.



vented cap

seal cap

Part no.	Description	packs/ctn
NKPETGF-0125S-A	conical shake flask 125ml, PETG material, seal cap, sterile	1pcs/bag, 24bag/ctn
NKPETGF-0125V-A	conical shake flask 125ml, PETG material, vented cap, sterile	1pcs/bag, 24bag/ctn
NKPCF-0125S-A	conical shake flask 125ml, PC material, seal cap, sterile	1pcs/bag, 24bag/ctn
NKPCF-0125V-A	conical shake flask 125ml, PC material, vented cap, sterile	1pcs/bag, 24bag/ctn
NKPETGF-0250S-A	conical shake flask 250ml, PETG material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPETGF-0250V-A	conical shake flask 250ml, PETG material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0250S-A	conical shake flask 250ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0250V-A	conical shake flask 250ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPETGF-0500S-A	conical shake flask 500ml, PETG material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPETGF-0500V-A	conical shake flask 500ml, PETG material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0500S-A	conical shake flask 500ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0500V-A	conical shake flask 500ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPETGF-1000V-A	conical shake flask 1000ml, PETG material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPETGF-1000S-A	conical shake flask 1000ml, PETG material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-1000S-A	conical shake flask 1000ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-1000V-A	conical shake flask 1000ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-3000V-A	big volume shake flask, 3000ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-5000V-A	big volume shake flask, 5000ml, PC material, vented cap, sterile	1pcs/bag, 6bag/ctn
NKPCF-3000S-A	big volume shake flask, 3000ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-5000S-A	big volume shake flask, 5000ml, PC material, seal cap, sterile	1pcs/bag, 6bag/ctn
NKKF-2800V-A	big volume shake flask, 2800ml, SBC material, vented cap, sterile	1pcs/bag, 6bag/ctn
NKKF-5000V-A	big volume shake flask, 5000ml, SBC material, vented cap, sterile	1pcs/bag, 4bag/ctn

Cell Culture Series Products



Part no.	Description	packs/ctn
NKPCF-0125S-B	baffled conical shake flask, 125ml, PC material, seal cap, sterile	1pcs/bag, 24bag/ctn
NKPCF-0125V-B	baffled conical shake flask, 125ml, PC material, vented cap, sterile	1pcs/bag, 24bag/ctn
NKPCF-0250S-B	baffled conical shake flask, 250ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0250V-B	baffled conical shake flask, 250ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0500S-B	baffled conical shake flask, 500ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-0500V-B	baffled conical shake flask, 500ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-1000S-B	baffled conical shake flask, 1000ml, PC material, seal cap, sterile	1pcs/bag, 12bag/ctn
NKPCF-1000V-B	baffled conical shake flask, 1000ml, PC material, vented cap, sterile	1pcs/bag, 12bag/ctn

Square medium bottle

Material Combination: Bottle body (PET/PETG/PC); Cap (HDPE/PP)

High Transparency: Glass-like clarity of the bottle body, with PETG version boasting 90% light transmittance

Practical Design: Precise and clear graduation marks; Square bottle shape for easy holding, storage and transportation

Safety & Protection: Crack-resistant, shatter-resistant and puncture resistant; Chamfered design of bottle mouth and cap ensures leak-proof performance after rigorous testing



60ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-60-A	PET square medium bottle, sterile	5x8, 5layer, 200pcs/ctn
NKSMB-60-B	PETG square medium bottle, sterile	5x8, 5layer, 200pcs/ctn
NKSMB-60-D-B	PETG square medium bottle, Ind wrap, sterile	1pcs/bag, 192pcs/ctn

30ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-30-A	PET square medium bottle, sterile	5x8, 7layer, 280pcs/ctn
NKSMB-30-B	PETG square medium bottle, sterile	5x8, 7layer, 280pcs/ctn

125ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-125-A	PET square medium bottle, sterile	4x6, 4layer, 96pcs/ctn
NKSMB-125-B	PETG square medium bottle, sterile	4x6, 4layer, 96pcs/ctn
NKSMB-125-D-B	PETG square medium bottle, Ind wrap, sterile	1pcs/bag, 80pcs/ctn
NKSMB-125-D-C	PC square medium bottle, Ind wrap, sterile	1pcs/bag, 200pcs/ctn

250ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-250-A	PET square medium bottle, sterile	5x6, 2layer, 60pcs/ctn
NKSMB-250-B	PETG square medium bottle, sterile	5x6, 2layer, 60pcs/ctn
NKSMB-250-D-B	PETG square medium bottle, Ind wrap, sterile	1pcs/bag, 60pcs/ctn
NKSMB-250-D-C	PC square medium bottle, Ind wrap, sterile	1pcs/bag, 100pcs/ctn

500ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-500-A	PET square medium bottle, sterile	4x5, 2layer, 40pcs/ctn
NKSMB-500-B	PETG square medium bottle, sterile	4x5, 2layer, 40pcs/ctn
NKSMB-500-D-B	PETG square medium bottle, Ind wrap, sterile	1pcs/bag, 48pcs/ctn
NKSMB-500-D-C	PC square medium bottle, Ind wrap, sterile	1pcs/bag, 100pcs/ctn


1000ml square medium bottle

Part no.	Description	packs/ctn
NKSMB-1000-A	PET square medium bottle, sterile	3x4, 2layer, 24pcs/ctn
NKSMB-1000-B	PETG square medium bottle, sterile	3x4, 2layer, 24pcs/ctn
NKSMB-1000-D-B	PETG square medium bottle, Ind wrap, sterile	1pcs/bag, 24pcs/ctn
NKSMB-1000-D-C	PC square medium bottle, Ind wrap, sterile	1pcs/bag, 50pcs/ctn

petri dish


Optically clear Petri Dishes made from Heavy Duty Food Grade Polystyrene plastic. Ethylene oxide sterilization, with lids, and have smooth surface and optimal transparent body. Optimized Design for Growing Bacteria - Perfectly sized 90mm dishes with triple vented lids allow an optimal airflow for cultivation of nutrient agar and to grow bacteria, mold, yeast, fungus and other microbes and organisms.

Part no.	Description	packs/ctn
NKPD-35	petri dish, 35mm, sterile	10pcs/bag, 200bags/ctn
NKPD-60	petri dish, 60mm, sterile	10pcs/bag, 100bags/ctn
NKPD-60-1	petri dish with grid, 60mm, sterile	10pcs/bag, 100bags/ctn
NKPD-90	petri dish, 90mm, sterile	10pcs/bag, 50bags/ctn
NKPD-90-2	petri dish, 90mm 2 room, sterile	10pcs/bag, 50bags/ctn
NKPD-90-3	petri dish, 90mm 3 room, sterile	10pcs/bag, 50bags/ctn
NKPD-90-4	petri dish, 90mm 4 room, sterile	10pcs/bag, 50bags/ctn
NKPD-100	petri dish, 100mm, sterile	10pcs/bag, 50bags/ctn
NKPD-100-1	petri dish, 100mm, square, sterile	10pcs/bag, 50bags/ctn
NKPD-120	petri dish, 120mm, sterile, 120*20mm	10pcs/bag, 50bags/ctn
NKPD-130	petri dish, 130mm, square, sterile, 130*17.5mm	300pcs/ctn
NKPD-150	petri dish, 150mm, sterile	5pcs/bag, 48bags/ctn



Inoculation loop/Cell Spreader

- Ultra-smooth loop surface design facilitates uniform and smooth streaking without damaging the gel surface
- Loops available in 1µL and 10µL sizes with strict measurement; the opposite end of each loop can be used as a needle.
- Round-shaped needle avoid any damage during colonies selection.



Part no.	Description	packs/ctn
NKIL-1-1	inoculation loop 1ul, Ind wrap, sterile	100pcs/pack, 5000pcs/ctn
NKIL-1-10	inoculation loop 1ul, 10pcs/bag, sterile	1000pcs/pack, 10000pcs/ctn
NKIL-1-20	inoculation loop 1ul, 20pcs/bag, sterile	2000pcs/pack, 10000pcs/ctn
NKIL-10-1	inoculation loop 10ul, Ind wrap, sterile	100pcs/pack, 5000pcs/ctn
NKIL-10-10	inoculation loop 10ul, 10pcs/bag, sterile	1000pcs/pack, 10000pcs/ctn
NKIL-10-20	inoculation loop 10ul, 20pcs/bag, sterile	2000pcs/pack, 10000pcs/ctn
NKS-L-1	spreader, L shape, Ind, sterile	200pcs/pack, 2000pcs/ctn
NKS-L-5	spreader, L shape, 5pcs/pack, sterile	500pcs/pack, 5000pcs/ctn
NKS-L-10	spreader, L shape, 10pcs/pack, sterile	500pcs/pack, 5000pcs/ctn
NKS-T-1	spreader, T shape, Ind, sterile	200pcs/pack, 2000pcs/ctn
NKS-T-5	spreader, T shape, 5pcs/pack, sterile	500pcs/pack, 5000pcs/ctn
NKS-T-10	spreader, T shape, 10pcs/pack, sterile	500pcs/pack, 5000pcs/ctn

Vacuum filter

- Vacuum-sealed and radiation-sterilized, DNase/RNase-free, non-pyrogenic, clean and safe.
- GL-45 threaded upper cup + universal hose connector, compatible with various containers and hose diameters.
- High transparency with clear graduations for easy volume visualization; filtered samples directly flow into sterile storage bottles.
- Optional membrane materials: NYLON/PVDF/PES; pore sizes: 0.22/0.45µm; multiple capacity options for upper cup/reservoir bottle.



Part no.	Description	packs/ctn
NKVF-150N-22	Vacuum filter 150ml, with 50mm 0.22um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-250N-22	Vacuum filter 250ml, with 50mm 0.22um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-500N-22	Vacuum filter 500ml, with 75mm 0.22um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-1000N-22	Vacuum filter 1000ml, with 91mm 0.22um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-150N-45	Vacuum filter 150ml, with 50mm 0.45um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-250N-45	Vacuum filter 250ml, with 50mm 0.45um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-500N-45	Vacuum filter 500ml, with 75mm 0.45um Nylon filter	1pcs/bag, 12bags/ctn
NKVF-1000N-45	Vacuum filter 1000ml, with 91mm 0.45um Nylon filter	1pcs/bag, 12bags/ctn



Part no.	Description	packs/ctn
NKVF-150P-22	Vacuum filter 150ml, with 50mm 0.22um PES filter	1pcs/bag, 12bags/ctn
NKVF-250P-22	Vacuum filter 250ml, with 50mm 0.22um PES filter	1pcs/bag, 12bags/ctn
NKVF-500P-22	Vacuum filter 500ml, with 75mm 0.22um PES filter	1pcs/bag, 12bags/ctn
NKVF-1000P-22	Vacuum filter 1000ml, with 91mm 0.22um PES filter	1pcs/bag, 12bags/ctn
NKVF-150P-45	Vacuum filter 150ml, with 50mm 0.45um PES filter	1pcs/bag, 12bags/ctn
NKVF-250P-45	Vacuum filter 250ml, with 50mm 0.45um PES filter	1pcs/bag, 12bags/ctn
NKVF-500P-45	Vacuum filter 500ml, with 75mm 0.45um PES filter	1pcs/bag, 12bags/ctn
NKVF-1000P-45	Vacuum filter 1000ml, with 91mm 0.45um PES filter	1pcs/bag, 12bags/ctn

Part no.	Description	packs/ctn
NKVF-150PV-22	Vacuum filter 150ml, with 50mm 0.22um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-250PV-22	Vacuum filter 250ml, with 50mm 0.22um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-500PV-22	Vacuum filter 500ml, with 75mm 0.22um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-1000PV-22	Vacuum filter 1000ml, with 91mm 0.22um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-150PV-45	Vacuum filter 150ml, with 50mm 0.45um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-250PV-45	Vacuum filter 250ml, with 50mm 0.45um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-500PV-45	Vacuum filter 500ml, with 75mm 0.45um PVDF filter	1pcs/bag, 12bags/ctn
NKVF-1000PV-45	Vacuum filter 1000ml, with 91mm 0.45um PVDF filter	1pcs/bag, 12bags/ctn



www.nokesci.com



WhatsAPP



Linkedin

JINHUA NOKE BIOTECHNOLOGY CO., LTD

Xianhua Technology Park, Qiubing St., Wucheng District, Jinhua, Zhejiang, China

Phone: +86 13373847308

Email: jane.tang@nokelab.com

